

DESCRIPTION OF
107 MICRODIALYSIS PUMP

The 107 Microdialysis Pump is specially designed for use with M Dialysis Catheters, 106 Syringe and Perfusion Fluid. The pump, syringe and catheter form an optimized system, where the microprocessor of the pump controls the high flush flow and the lower normal flow. The normal flow is adjustable in eight different flow rates.

Intended use: The 107 Microdialysis Pump is a portable syringe pump intended to pump perfusion fluid and dialysate through a microdialysis catheter, enabling microdialysis sampling.

The pump is intended to be operated by medically trained staff.

Note! The 107 Microdialysis Pump must only be used for microdialysis, together with microdialysis catheters and accessories from M Dialysis AB.





Fig. 1. A: Pump: (1) Cover (2) Battery Cover.

DESCRIPTION OF
107 MICRODIALYSIS PUMP

Fig 1. B: Pump from above: (3) Light-emitting diodes (4) Power switch (5) Drive screw (6) Flow Change Switch (7) Tool for changing the flow

Fig 1. C: Syringe: (8) Piston with thread



BATTERY CHANGE




Fig. 2. A: Battery cover closed. B: Battery cover open. C: Battery with pull tab

1. Open the battery cover by pressing with your thumb in the direction of the arrow (see Fig. 2 A,B).
2. Remove the old battery by pulling the label tab.
3. Attention: Be sure to insert the battery properly. Insert the positive end (+) of the battery to the positive end of plate first, followed by the negative (-) side to avoid a short circuit.
4. Fold the label down.
5. Replace the battery cover.

Low battery is indicated by two red light signals every ten seconds. Battery replacement can be made during a microdialysis investigation which is underway provided that the syringe remains in the pump and that the cover is not opened.

Used batteries should be disposed of according to local environmental regulations or contact M Dialysis AB for more information.

Note! If the pump is not in use, the batteries should be removed.

CONNECTION OF SYRINGE




Fig. 3. Positioning of syringe in pump.

1. Fill the syringe with a maximum of 2.5 mL Perfusion Fluid at room temperature
2. Remove any air bubbles. Remove the filling needle from the syringe.
3. Connect the luer connector of the syringe to the inlet tubing of the microdialysis catheter. Before flushing, fit a microvial into the microvial holder.
4. Fit the syringe into the pump by inserting the front section first (1). Allow the syringe to drop into position so that the threaded piston (2) lines up with the drive screw (see Fig 3).
5. Close the cover by pressing on the sides of the cover so there is no "click". The pump will start automatically and the green light will flash every other second during the flush sequence (5 min). Note that each time the pump is opened and closed with a syringe in place; a new flushsequence will start.
6. Check that there is fluid in the microvial after the flush sequence and change to a new vial. If there is no fluid in the microvial, start a new flush sequence by opening and closing the pump cover.

PUMP FUNCTIONS

Standby: No syringe in pump.

Flush sequence: Syringe in the pump. High flow for five minutes (15µL/min), normal flow thereafter.

Normal flow: Adjustable in eight flow rates of 0, 0.1, 0.2, 0.3, 0.5, 1, 2, 5 µL/min.

To indicate the current function of the pump, there are two light-emitting diodes (LEDs, see fig. 1) located under the transparent cover by the "nose" of the pump. The LEDs can illuminate red or green.

When the syringe is inserted into the pump and the cover is closed, a switch underneath the front end of the syringe is activated (see Fig. 1) and the flush sequence starts.

During the flush sequence, a green light will flash every other second (see Light signals). The flush sequence is programmed so that the complete microdialysis catheter is filled with fluid and all air bubbles are removed. When the flush is complete the pump changes to normal flow which is indicated with a green light flash every 10 seconds (see Light signals).


ADJUSTING NORMAL FLOW

The normal flow can be adjusted in eight different flow rates by the switch at the back of the pump.

1. Take out the tool from the bottom of the pump. (See Fig. 4 A.)
2. Put the tool in the arrow of the flow rate switch. (See Fig. 4 B.)
3. Turn the switch until the arrow lines up with the desired flow rate.
4. After a few seconds a light signal confirms that a flow rate change has occurred. (See Light signals.)

The flow can be adjusted during a running microdialysis investigation if needed.

Fig. 4. A,B



LIGHT SIGNALS

INFORMATION SIGNALS

System OK
After correct insertion of battery: Three green/red light signals ending in one green light signal of about 3 seconds.

Flush sequence
Green light signal every other second.

Normal flow
Green light signal every 10 seconds.

Flow change
Three green light signals with variation of duration depending on the chosen normal flow (see accompanying 106/107 Pump Safety Test Sheet - 8050067).

ERROR SIGNALS

Pump error
Red light signal every 5 seconds if the problem is during flushing, every 10 seconds if the problem is during normal flow.

Low battery
Two red light signals every 10 seconds.

The error signal stops when the fault has been corrected (e.g. filling the syringe, removing the blockage in the syringe, changing the battery etc.)

CLEANING

Use a moist soft cloth to clean the casing of the pump.
Do not use abrasive cleaners.
Recommended cleaning substances:
- Soap solution
- 70 % ethanol

WATER RESISTANCE

The pump is splash proof and can tolerate short-duration splashing with water.
This means that the patient can shower with the pump functioning if the pump is protected with a plastic bag. If the pump is accidentally immersed in water, the syringe and battery compartments must be dried out.

WARRANTY

M Dialysis AB guarantees all components of the 107 Microdialysis Pump to be free from defects of material and workmanship for a period of 12 months after initial purchase.
M Dialysis AB will repair or replace, at our discretion, the 107 Microdialysis Pump during the aforementioned warranty period.
M Dialysis AB reserves the right to waive all warranties in the case of problems due to improper handling, improper field of application or unauthorised modifications.
For warranty repair, the 107 Microdialysis Pump must be returned to M Dialysis AB or to an authorised representative. The Owner shall prepay shipping charges to M Dialysis, and M Dialysis AB shall pay shipping charges to return the product to the Owner.

WARNINGS

- No modification or repair of this equipment is allowed.
- Do not drop this equipment, consult the safety and control instruction.
- Only use battery from M Dialysis, ref no 8001788
- Connect only items that have been specified as part of the equipment.
- Remove battery from equipment when not in use
- Never submerge the equipment in water or other liquid
- Battery compartment shall be clean and dry.
- Read instruction for battery handling on battery package labeling
- For battery disposal follow hospital procedure for disposing of batteries.

- In order not to affect the function of the pump, ensure that radiotransmitters, mobile telephones and other wireless communication equipment is used at a safe distance from the 107 Microdialysis Pump.
- 107 Microdialysis Pump should not be exposed to disturbance levels exceeding those given in IEC 60601-1-2.
- If the 107 Microdialysis Pump is to be discarded, please contact M Dialysis AB or your local supplier for more information.
- Remove the pump before MRI scanning.
- To avoid skin irritation, do not place the pump in direct contact with the skin.

EXPLANATION OF SYMBOLS

IPX4

Protected against splashing water.



Follow instructions for use.



Pump fulfills requirements of IEC 601, safety class CF.

CE 0413

Certified according to the Medical Device Directive, (Intertek, Sweden)




Separate collection for waste of electrical and electronic equipment.




ETL - listed product. Conforms to UL 60601-1:2003 Rev 2006 and CSA C22.2#601.1 (R2001)



MR Unsafe



Manufacturer and date of manufacture



REF

Catalogue number

SN

Serial number


CHECK FLOW SETTING!

Check the Normal Flow setting before use.

TECHNICAL SPECIFICATION

Dimension:	90 x 50 x 20 mm
Weight:	70 g (incl. battery)
Battery:	6V Silver oxide Ref 8001788
Battery lifetime:	10 days at flow of 0.3 µL/min
Casing:	ABS plastic, splash-proof
Flush flow:	15 µL/min
Normal flow:	0-5 µL/min
Operating temperature:	+5 - +40 °C
Fehlersignal:	Pump error Low battery CE
107	CE
106 Syringe	CE
Oxygen Rich Environment:	The pump is not tested in a Oxygen Rich Environment.

107 MICRODIALYSIS PUMP



CONTENTS

Description of Pump	2,3
Battery Change	4
Connection of the Syringe	5
Pump Functions	6
Adjusting Normal Flow	7
Light Signals	8
Cleaning, Water Resistance and Warranty	9
Warnings	10
Explanation of Symbols	11
Technical Specification	12

μdialysis

Manufactured by:
M Dialysis AB
Hammarby Fabriksv 43 SE-120 30 • Stockholm • Sweden
Tel: +46-8-470 10 20
E-mail: info@mdialysis.com

USA office:
73 Princeton Street
N.Chelmsford • MA 01863 • USA
Phone: +1 978 251 1940, +1-866-868-9236
Fax: +1 978 251-1960
E-mail: usa@mdialysis.com