REAGENT Set A

For ISCUS/ISCUS^{flex} Clinical Microdialysis Analyzer

Ref. No. 8002163

Content

1. Reagent: One bottle lyophilized reagent each for Glucose, Lactate, Pyruvate and Glycerol placed in a cassette.
2. Buffer: One bottle 6 mL each for Glucose, Lactate, Pyruvate and Glycerol.
3. Calibrator: One bottle of Calibrator A 6 mL placed in the Reagent Cassette.

Preparation

1. Unscrew the cap with the membrane from the reagent and calibrator bottles, placed in the cassette. Remove and discard the rubber stoppers.
2. Unscrew the cap from the buffer bottles and gently transfer the contents to the corresponding reagent bottle.
3. Fasten the cap with the membrane on the reagent and calibrator bottles in the cassette, without returning the rubber stoppers.
4. Let the reagents and calibrator stand and equilibrate in room temperature for at least 30 minutes prior to use. The contents will be completely mixed when the Reagent Cassette is placed and identified in the analyzer.

Intended use and stability of solution


The Reagent Set is made for use in ISCUS or ISCUS^{flex} Clinical Microdialysis Analyzer. Each Reagent Cassette has a unique code, written on the label, which should be registered when placing it in the analyzer.

The Reagent Set is stable up to expiry date when stored at +2 to +8 °C. Reconstituted reagents are stable for five days in the analyzer. The contents are sufficient for around 310 analyses.

GLUCOSE

Colorimetric method for the quantitative determination of Glucose in Microdialysates.

Measuring principle
Glucose is enzymatically oxidized by glucose oxidase (GOD). The hydrogen peroxide formed reacts with phenol and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields a red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the glucose concentration.

Linear range: 0.1 - 25 mmol/L

	Component in test solution	Concentration
Glucose reagent	4-aminoantipyrine	0.77 mmol/L
	Ascorbate oxidase	>3 kU/L
	Glucose oxidase	>1.5 kU/L
	Peroxidase	>1.5 kU/L
Glucose buffer	Phosphate buffer, pH 7.0	0.1 mol/L
	Phenol	11 mmol/L
	Sodium azide	0.4 g/L



References:

1. Barhem and P. Trinder, Analyst 97 (1972) 142

LACTATE

Colorimetric method for the quantitative determination of Lactate in Microdialysates.

Measuring principle
Lactate is enzymatically oxidized by lactate oxidase. The hydrogen peroxide formed reacts with 4-chlorophenol and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields a red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the lactate concentration.

Linear range: 0.1 - 12 mmol/L

	Component in test solution	Concentration
Lactate reagent	4-aminoantipyrine	0.4 mmol/L
	Lactate oxidase	>500 U/L
	Peroxidase	>500 U/L
	Ascorbate oxidase	>12.0 kU/L
Lactate buffer	PIPES buffer, pH 6.8	1.00 mmol/L
	4-Chlorophenol	5.4 mmol/L
	Sodium oxalate	7.5 mmol/L
	EDTA-disodium salt	5 mmol/L
	Sodium azide	0.3 g/L


References:

1. N. Shimajo et al. Clin Chem 35 (1989) 1992
2. H.F. Kühnle et al. J.Clin Chem BioChem 15 (1977) 171
3. T.O. Kleine et al. Dtsch Med Wschr 104 (1979) 553

PYRUVATE

Colorimetric method for the quantitative determination of Pyruvate in Microdialysates.

Measuring principle
Pyruvate is enzymatically oxidized by pyruvate oxidase (PyrOx). The hydrogen peroxide formed reacts with N-ethyl-N-(2-hydroxy-3-sulfopropyl)-m-toluidine and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields a red-violet colored quinonediimine. The rate of formation is measured photometrically at 530 nm and is proportional to the pyruvate concentration.

Linear range: 2 - 300 (10 - 1500) µmol/L

	Component in test solution	Concentration
Pyruvate reagent	4-amino-antipyrine	0.3 mmol/L
	Tiamine pyrophosphate	0.2 mmol/L
	FAD	10 µmol/L
	Pyruvate Oxidase	>0.2 kU/L
	Peroxidase	>0.8 kU/L
	Ascorbate Oxidase	>10 kU/L
Pyruvate buffer	Citrate buffer, pH 6.1	100 mmol/L
	Potassium dihydrogenphosphate	10 mmol/L
	MgCl2	10 mmol/L
	TOOS	1.5 mmol/L
	Sodium azide	0.3 g/L

References:



1. B. Sedewitz, et al., J. Bacteriol, 160 (1984) 273-278
2. M. Nawata, et al., Anal Biochem, 190 (1990) 84-87
3. H. Araki and M. Yamada, in: H. U. Bergmeyer (Editor), Methods of Enzymatic Analysis, 3rd ed., Vol 6, Verlag Chemie, Weinheim, 1984

GLYCEROL

Colorimetric method for the quantitative determination of Glycerol in Microdialysates.

Measuring principle

Glycerol is phosphorylated by adenosin triphosphate (ATP) and glycerol kinase (GK) to glycerol-3-phosphate, which is subsequently oxidized in the presence of glycerol-3-phosphate oxidase (GPO). The hydrogen peroxide formed reacts with 3,5-dichloro-2-hydroxy-benzene sulphonic acid (DCHBS) and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields a red-violet colored quinoneimine (ACBS). The rate of formation is measured photometrically at 530 nm and is proportional to the glycerol concentration.

Linear range: 10 - 1500 µmol/L

	Component in test solution	Concentration
Glycerol reagent	4-aminoantipyrine	0.4 mmol/L
	ATP	1.0 mmol/L
	Glycerol kinase	>400 U/L
	Glycerol-3-phosphate-oxidase	>1.5 kU/L
	Peroxidase	>1 kU/L
Glycerol buffer	Ascorbate oxidase	>7.0 kU/L
	PIPES buffer, pH 7.6	40 mmol/L
	DCHBS	1.5 mmol/L
	Magnesium ions	17.5 mmol/L
	Sodium azide	0.2 g/L

References:

1. K.J. Foster and K.G.M.M. Alberti, Clin Chem 24 (1978) 1568

CALIBRATOR A

Calibration values	
Glucose	5.55 mmol/L
Lactate	2.5 mmol/L
Glycerol	475 µmol/L
Pyruvate	250 µmol/L

Sample material	WARNING
Microdialysates	Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.
For in vitro use only	The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
Symbol declaration	Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10 % sodium hydroxide.
 Last day of use	
 Lot number	
 Storage temperature	
 The product meets EU directive for IVD (98/79/EC)	

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M Dialysis AB