REAGENT Set B

For ISCUS/ISCUS^{flex} Clinical Microdialysis Analyzer

Ref. No. 8002164

Content

1. Reagent: One bottle lyophilized reagent each for Glucose, Lactate, Pyruvate, Glycerol and Glutamate placed in a cassette.
2. Buffer: One bottle 6 mL each for Glucose, Lactate, Pyruvate and Glycerol and 4 mL for Glutamate.
3. Calibrator: One bottle of Calibrator A 6 mL placed in the Reagent Cassette.

Preparation

1. Unscrew the cap with the membrane from the reagent and calibrator bottles, placed in the cassette. Remove and discard the rubber stoppers.
2. Unscrew the cap from the buffer bottles and gently transfer the contents to the corresponding reagent bottle.
3. Fasten the cap with the membrane on the reagent and calibrator bottles in the cassette, without returning the rubber stoppers.
4. Let the reagents and calibrator stand and equilibrate in room temperature for at least 30 minutes prior to use. The contents will be completely mixed when the Reagent Cassette is placed and identified in the analyzer.

Intended use and stability of solution


The Reagent Set is made for use in ISCUS or ISCUS^{flex} Clinical Microdialysis Analyzer. Each Reagent Cassette has a unique code, written on the label, which should be registered when placing it in the analyzer.

The Reagent Set is stable up to expiry date when stored at +2 to +8 °C. Reconstituted reagents are stable for five days in the analyzer. The contents are sufficient for around 215 analyses.

GLUCOSE

Colorimetric method for the quantitative determination of Glucose in Microdialysates.

Measuring principle
Glucose is enzymatically oxidized by glucose oxidase (GOD). The hydrogen peroxide formed reacts with phenol and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields a red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the glucose concentration.

Linear range: 0.1 - 25 mmol/L

	Component in test solution	Concentration
Glucose reagent	4-aminoantipyrine Ascorbate oxidase Glucose oxidase Peroxidase	0.77 mmol/L >3 kU/L >1.5 kU/L >1.5 kU/L
Glucose buffer	Phosphate buffer, pH 7.0 Phenol Sodium azide	0.1 mol/L 11 mmol/L 0.4 g/L



References:

1. Barhem and P. Trinder, Analyst 97 (1972) 142

LACTATE

Colorimetric method for the quantitative determination of Lactate in Microdialysates.

Measuring principle
Lactate is enzymatically oxidized by lactate oxidase. The hydrogen peroxide formed reacts with 4-chlorophenol and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields a red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the lactate concentration.

Linear range: 0.1 - 12 mmol/L

	Component in test solution	Concentration
Lactate reagent	4-aminoantipyrine Lactate oxidase Peroxidase	0.4 mmol/L >500 U/L >500 U/L
Lactate buffer	Ascorbate oxidase PIPES buffer, pH 6.8 4-Chlorophenol Sodium oxalate EDTA-disodium salt Sodium azide	>12.0 kU/L 100 mmol/L 5.4 mmol/L 7.5 mmol/L 5 mmol/L 0.3 g/L



References:

1. N. Shimojo et al. Clin Chem 35 (1989) 1992
2. H.F. Kühnle et al. J.Clin Chem BioChem 15 (1977) 171
3. T.O. Kleine et al. Dtsch Med Wschr 104 (1979) 553

PYRUVATE

Colorimetric method for the quantitative determination of Pyruvate in Microdialysates.

Measuring principle
Pyruvate is enzymatically oxidized by pyruvate oxidase (PyrOx). The hydrogen peroxide formed reacts with N-ethyl-N-(2-hydroxy-3-sulfo-propyl)-m-toluidine and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields a red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the pyruvate concentration.

Linear range: 2 - 300 (10 - 1500) µmol/L

	Component in test solution	Concentration
Pyruvate reagent	4-amino-antipyrine Thiamine pyrophosphate FAD Pyruvate Oxidase Peroxidase	0.3 mmol/L 0.2 mmol/L 10 µmol/L >0.2 kU/L >0.8 kU/L
Pyruvate buffer	Ascorbate Oxidase Citrate buffer, pH 6.1 Potassium dihydrogenphosphate MgCl ₂ TOOS Sodium azide	>10 kU/L 100 mmol/L 10 mmol/L 10 mmol/L 1.5 mmol/L 0.3 g/L

References:



1. B. Sedewitz, et al., J. Bacteriol, 160 (1984) 273-278
2. M. Nawata, et al., Anal Biochem, 190 (1990) 84-87
3. H. Araki and M. Yamada, in: H. U. Bergmeyer (Editor), Methods of Enzymatic Analysis, 3rd ed., Vol 6, Verlag Chemie, Weinheim, 1984

GLYCEROL

Colorimetric method for the quantitative determination of Glycerol in Microdialysates.

Measuring principle

Glycerol is phosphorylated by adenosin triphosphate (ATP) and glycerol kinase (GK) to glycerol-3-phosphate, which is subsequently oxidized in the presence of glycerol-3-phosphate oxidase (GPO). The hydrogen peroxide formed reacts with 3,5-dichloro-2-hydroxy-benzene sulphonic acid (DCHBS) and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields a red-violet colored quinoneimine (ACBS). The rate of formation is measured photometrically at 530 nm and is proportional to the glycerol concentration.

Linear range: 10 - 1500 µmol/L

	Component in test solution	Concentration
Glycerol reagent	4-aminoantipyrine ATP Glycerol kinase Glycerol-3-phosphate-oxidase Peroxidase	0.4 mmol/L 1.0 mmol/L >400 U/L >1.5 kU/L >1 kU/L
Glycerol buffer	Ascorbate oxidase PIPES buffer, pH 7.6 DCHBS Magnesium ions Sodium azide	>7.0 kU/L 40 mmol/L 1.5 mmol/L 17.5 mmol/L 0.2 g/L

References:



1. K.J. Foster and K.G.M.M. Alberti, Clin Chem 24 (1978) 1568

GLUTAMATE

Colorimetric method for the quantitative determination of Glutamate in Microdialysates.

Measuring principle

Glutamate is enzymatically oxidized by glutamate oxidase (GtOx). The hydrogen peroxide formed reacts with N-ethyl-N-(2-hydroxy-3-sulfo-propyl)-m-toluidine (TOOS) and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields the red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the glutamate.

Linear range: 1 - 150 µmol/L

	Component in test solution	Concentration
Glutamate reagent	4-amino-antipyrine Glutamate Oxidase Peroxidase Ascorbate Oxidase	0.3 mmol/L >0.25 kU/L >0.8 kU/L >17.5 kU/L
Glutamate buffer	PIPES buffer, pH 6.8 TOOS Sodium azide	100 mmol/L 1.5 mmol/L 0.3 g/L

References:

1. H. Kusakabe et al., Agric. Biol. Chem., 47 (1983) 179
2. H. Kusakabe et al., Agric. Biol. Chem., 47 (1983) 1323
3. A. Böhrner et al., Eur. J. Biochem. 182 (1989) 327-332
4. U. Wollenberger et al., Biosensors 4 (1989) 381-39

CALIBRATOR A

Calibration values

Glucose	5.55 mmol/L
Lactate	2.5 mmol/L
Glycerol	475 µmol/L
Pyruvate	250 µmol/L
Glutamate	25 µmol/L

Sample material	Microdialysates	WARNING
Symbol declaration	 Last day of use	Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents. The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention. Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10 % sodium hydroxide.
 LOT	Lot number	
	Storage temperature	
	Consult instructions for use	
	In vitro diagnostic device	For in vitro use only
	The product meets EU directive for IVD (98/79/EC)	

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M Dialysis AB