REAGENT KIT

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Referenser:
1. Barhem and P.Trinder, Analyst 97(1972)142



REAGENT KIT

for the 600 and ISCUS^{flex} Microdialysis Analyzers

GLUCOSE

Colorimetric method for the quantitative determination of Glucose in Microdialysates.

Measuring principle
Glucose is enzymatically oxidised by glucose oxidase (GOD). The hydrogen peroxide formed reacts with phenol and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields the red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the glucose concentration.

Linear range: 0.1 - 25 mmol/L

	Component	Concentration in test solution
Glucose reagent	4-aminoantipyrine Ascorbate oxidase Glucose oxidase	0.77 mmol/L >3 kU/L >1.5 kU/L
Glucose buffer	Peroxidase Phosphate buffer, pH 7.0 Phenol Sodium azide	>1.5 kU/L 0.1 mol/L 11 mmol/L 0.4 g/L

Sample material Microdialysates	WARNING: Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.
For in vitro use only	The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
Symbol declaration:	Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10 % sodium hydroxide
 Last day of use	
 Lot number	
 Storage temperature	
 The product meets EU directive for IVD (98/79/EC)	

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

IMPORTANT!
Remove the rubber stopper under the membrane lid before placing the calibrator in the Microdialysis Analyzer.

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



for the 600 and ISCUS^{flex} Microdialysis Analyzers

GLYCEROL

Colorimetric method for the quantitative determination of Glycerol in Microdialysates.

Measuring principle
Glycerol is phosphorylated by adenosine triphosphate (ATP) and glycerol kinase (GK) to glycerol-3-phosphate, which is subsequently oxidized in the presence of glycerol-3-phosphate oxidase (GPO). The hydrogen peroxide formed reacts with 3,5-dichloro-2-hydroxy-benzene sulphonic acid (DCHBS) and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields the red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the glycerol concentration.

Linear range: 0.01 - 1.5 mmol/L

	Component	Concentration in test solution
Glycerol reagent	4-aminoantipyrine ATP Glycerol kinase Glycerol-3-phosphate-oxidase Peroxidase Ascorbate oxidase	0.4 mmol/L 1.0 mmol/L >400 U/L >1.5 kU/L >1 kU/L >7.0 kU/L
Glycerol buffer	PIPES buffer, pH 7.6 DCHBS Magnesium ions Sodium azide	40 mmol/L 1.5 mmol/L 17.5 mmol/L 0.2 g/L

Sample material Microdialysates	WARNING: Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.
For in vitro use only	The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
Symbol declaration:	Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10 % sodium hydroxide
 Last day of use	
 Lot number	
 Storage temperature	
 The product meets EU directive for IVD (98/79/EC)	

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


Preparation and stability of solution

- Unscrew the cap with the membrane from the reagent bottle. Remove and discard the rubber stopper.
- Transfer the contents of the buffer bottle to the reagent bottle.
- Fasten the cap with the membrane on the reagent bottle, without Rubber stopper.
- Dissolve contents completely by gently turning the bottle upside-down at least ten times. Let the reagent stand and equilibrate in room temperature for at least 30 minutes prior to use. Reconstituted reagent is stable for five days in the instrument.

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

for the 600 and ISCUS^{flex} Microdialysis Analyzers

LACTATE

Colorimetric method for the quantitative determination of Lactate in Microdialysates.

Measuring principle
Lactate is enzymatically oxidised by lactate oxidase. The hydrogen peroxide formed reacts with 4-chlorophenol and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields the red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the lactate concentration.

Linear range: 0.1 - 12 mmol/L

	Component	Concentration in test solution
Lactate reagent	4-aminoantipyrine Lactate oxidase Peroxidase Ascorbate oxidase	0.4 mmol/L >500 U/L >500 U/L >12.0 kU/L
Lactate buffer	PIPES buffer, pH 6.8 4-Chlorophenol Sodium oxalate EDTA-disodium salt Sodium azide	100 mmol/L 5.4 mmol/L 7.5 mmol/L 5 mmol/L 0.3 g/L

Sample material Microdialysates	WARNING: Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.
For in vitro use only	The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
Symbol declaration:	Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10 % sodium hydroxide
 Last day of use	
 Lot number	
 Storage temperature	
 The product meets EU directive for IVD (98/79/EC)	

References:
1.N.Shimojo et al, Clin Chem 35(1989)1992 2.H.F.Kühnle et al, J.Clin Chem BioChem 15 (1977)171 3.T.O.Kleine et al, Dtsch Med Wscht 104 (1979) 553

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
GB

- Reagent: One bottle lyophilized reagent each for glucose, lactate, pyruvate and glycerol
- Buffer: One bottle à 6 mL each for glucose, lactate, pyruvate and glycerol
- Calibrator: One bottle à 6 mL
Reagents are sufficient for 350 determinations.
Reagents and calibrator are stable up to expiry date when stored at +2 to +8 °C

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

for the 600 and ISCUS^{flex} Microdialysis Analyzers

PYRUVATE

Colorimetric method for the quantitative determination of Pyruvate in Microdialysates.

Measuring principle
Pyruvate is enzymatically oxidized by pyruvate oxidase (PyrOx). The hydrogen peroxide formed reacts with N-ethyl-N-(2-hydroxy-3-sulfopropyl)-m-toluidine and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields the red-violet colored quinonedimine. The rate of formation is measured photometrically at 530 nm and is proportional to the pyruvate concentration.

Linear range: 2 - 300 (10 - 1500) µmol/L

	Component	Concentration in test solution
Pyruvate reagent	4-amino-antipyrine Tiamine pyrophosphate FAD Pyruvate Oxidase Peroxidase Ascorbate Oxidase	0.3 mmol/L 0.2 mmol/L 10 µmol/L >0.2 kU/L >0.8 kU/L >10 kU/L
Pyruvate buffer	Citrate buffer, pH 6.1 Potassium dihydrogenphosphate MgCl ₂ TOOS Sodium azide	100 mmol/L 10 mmol/L 10 mmol/L 1.5 mmol/L 0.3 g/L

Sample material Microdialysates	WARNING: Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.
For in vitro use only	The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
Symbol declaration:	Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10 % sodium hydroxide
 Last day of use	
 Lot number	
 Storage temperature	
 The product meets EU directive for IVD (98/79/EC)	

References:
1. B. Sedewitz, et al., J. Bacteriol., 160 (1984) 273-278
2. M. Nawata, et al., Anal Biohem., 190 (1990) 84-87
3. H. Araki and M. Yamada, in: H. U. Bergmeyer (Editor), Methods of Enzymatic Analysis, 3rd ed., Vol 6, Verlag Chemie, Weinheim, 1984

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Preparation and stability of solution

- Unscrew the cap with the membrane from the reagent bottle. Remove and discard the rubber stopper.
- Transfer the contents of the buffer bottle to the reagent bottle.
- Fasten the cap with the membrane on the reagent bottle, without Rubber stopper.
- Dissolve contents completely by gently turning the bottle upside-down at least ten times. Let the reagent stand and equilibrate in room temperature for at least 30 minutes prior to use.
Reconstituted reagent is stable for five days in the instrument.

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