REAGENS

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UREA

REF: P000026


UV-method for the quantitative determination of urea in Microdialysates.

Also required:

Calibrator A. REF: P000057

Measuring principle

Urea is hydrolyzed in the presence of urease to ammonium ions and carbon dioxide. The ammonium ions react with 2-oxoglutarate in the presence of glutamate dehydrogenase (GIDH) and NADH to form glutamate and NAD⁺. The rate of utilization of NADH is measured photometrically at 375 nm and is proportional to the urea concentration.

Linear range: 0.5 - 25 mmol/L (recently prepared), 0.5 - 17 mmol/L (after three days)

Reagents

1. Reagent: 5 bottles of lyophilisate
2. Buffer: 5 bottles of 6 mL.

Reagent sufficient for 5 x 350 determinations.

Reagents are stable up to expiry date when stored at +2 to +8°C.

Preparation and stability of solution

1. Unscrew the cap with the membrane from the reagent bottle. Remove and discard the rubber stopper.
2. Transfer the contents of the buffer bottle to the reagent bottle.
3. Fasten the cap with the membrane on the reagent bottle, without Rubber stopper.
4. Dissolve contents completely by gently turning the bottle upside-down at least ten times. Let the reagent stand and equilibrate in room temperature for at least 30 minutes prior to use.

Reconstituted reagent is stable for three days in the instrument.

Component	Concentration in test solution
Urea reagent	Urease
	Glutamate dehydrogenase
	NADH
	2-oxoglutarate
Urea buffer	TRIS buffer, pH 8.0
	Sodium azide
	0.3 g/L

Sample material Microdialysates	Calibration Use Calibrator A. REF: P000057
------------------------------------	---

For in vitro use only. **WARNING:**

Symbol declaration:

 Last day of use

 Lot number

 Storage temperature

 See instructions for use

 In vitro diagnostic reagent

 The product meets EU directive for IVD (98/79/EC)

Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.
The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10 % sodium hydroxide

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
Quality Control:

It is recommended that the system be controlled using Control Samples. These controls should be assayed according to local quality assurance schemes. The use of other control materials has not been evaluated. Satisfactory level of performance is achieved when the analyte values for the controls are within the "Acceptable Control Range" published in the Packet Insert with the controls.

Additional information is found in the Technical Manual for the Microdialysis Analyzers.

Manufactured by:
M Dialysis AB
Hammarby Fabriksväg 43
SE-120 30 • Stockholm • Sweden
Tel: +46-8-470 10 20
Fax: +46-8-470 10 55
E-mail: info@mdialysis.com

USA office:
73 Princeton Street
N.Chelmsford • MA 01863 • USA
Phone: +1 978 251 1940, +1-866-868-9236
Fax: +1 978 251-1960
E-mail: usa@mdialysis.com