


## Early detection of secondary insults & software for Multimodal Monitoring

The Consensus paper (Hutchinson et al. Intensive Care Med. 2015) from the International Microdialysis Forum meeting in Cambridge includes; experts guidance for use of microdialysis in traumatic brain injury and subarachnoid hemorrhage, catheter location, reference values and interventions. Find below a few conclusions from the paper:

- Low brain glucose is associated with unfavorable outcome.
- An increase in the Lactate Pyruvate (LP) ratio in the presence of low pyruvate (and low oxygen) indicates ischemia.
- An increase in the LP ratio in the presence of normal or high pyruvate (and normal oxygen) indicates mitochondrial dysfunction.
- A high LP ratio is associated with unfavorable outcome.

### ICUpilot - software for multimodal monitoring

ICUpilot is a unique tool for multimodal monitoring in the ICU. Bedside Patient Monitors (showing e.g. pulse, blood pressure, ICP, CPP) as well as the Microdialysis Analyzer can be connected to a separate computer for on-line analysis and comparison of all data collected bedside during the entire care of the patient.



## Literature:

**Consensus statement from the 2014 International Microdialysis Forum.** Hutchinson PJ et al. Intensive Care Med. 2015 Sep;41(9):1517-28.

**Implementation of cerebral microdialysis at a community-based hospital: A 5-year retrospective analysis.** Surg Neurol Int. 2012  
Chen et al. Department of Neurosurgery, Legacy Emanuel Medical Center, Portland, USA.

**International Multidisciplinary Consensus Conference on Multimodality Monitoring: Cerebral Metabolism.**  
Hutchinson P, O'Phelan K; The Participants in the International Multidisciplinary Consensus Conference on Multimodality Monitoring. Neurocrit Care. 2014 Sep 11

**Intracerebral monitoring of silent infarcts after subarachnoid hemorrhage.** Neurocrit Care. 2011 Apr;14(2):162-7.  
Helbok et al. Division of Critical Care Neurology, Department of Neurology, Columbia University Medical Center, Milstein Hospital 8 Center 177 Fort Washington Ave, New York, NY, USA

**Neuromonitoring with microdialysis in severe traumatic brain injury patients.** Sanchez JJ, Bidot CJ, O'Phelan K, Gajavelli S, Yokobori S, Olvey S, Jagid J, Garcia JA, Nemeth Z, Bullock R. Acta Neurochir Suppl. 2013  
University of Miami Miller School of Medicine, Miami, FL 33136, USA

**Neuromonitoring in Intensive Care: Focus on Microdialysis and Its Nursing Implications,** 2009 American Association of Neuroscience Nurses  
Mary Presciutti, J. Michael Schmidt, Sheila Alexander

**Cerebral Metabolism and the Role of Glucose Control in Acute Traumatic Brain Injury.** Buitrago Blanco MM, Prashant GN, Vespa PM. Neurosurg Clin N Am. 2016 Oct;27

### M Dialysis AB

M Dialysis is the leading company devoted to the development, manufacturing and marketing of the Microdialysis technique.  
The head office is located in Stockholm, Sweden, with a subsidiary in Boston MA, USA. M Dialysis has distributors across the globe, responsible for local sales, service and support.

**mdialysis**

M Dialysis AB, Box 5049, SE-121 05 Stockholm, Sweden,  
Tel: +46 8 470 10 20, Fax: +46-8-470 10 55  
E-mail: info@mdialysis.se

73 Princeton Street, North Chelmsford, MA 01863, USA  
Tel: (866) 868-9236, (978) 251-1940 Fax: (978) 251-1960  
E-mail: usa@mdialysis.com  
[www.mdialysis.com](http://www.mdialysis.com)



# Microdialysis

## Advanced monitoring in neurointensive care




Ref. No. 8011009F 2017 08 17

**mdialysis**

## Microdialysis in Neurointensive Care

Microdialysis is a tool for *in vivo* sampling of soft tissues that utilizes the principal of diffusion through a semi-permeable membrane. The technology is minimally-invasive, easy to handle, and may be used continuously over a period of several days.

The method is performed by inserting a Microdialysis catheter into the tissue being studied. The Microdialysis membrane of the catheter is in direct contact with the soft tissue.



The catheter is perfused with a sterile isotonic solution via a small pump attached to its inlet lumen. In the tissue, substances from the interstitial fluid diffuse through the semi-permeable Microdialysis membrane into the perfusion fluid. This fluid, now known as dialysate, moves through the outlet lumen and into a collection microvial. Microvials are exchanged at regular intervals. The dialysate collected may be analyzed immediately using the ISCUS<sup>flex</sup> Microdialysis Analyzer.


The metabolite values in the dialysate provide a picture of the local tissue metabolism. This has been particularly useful in neurointensive care as there are well described metabolic changes that occur with secondary ischemic events in the cases of traumatic brain injury (TBI) and subarachnoid hemorrhage (SAH).

**Secondary ischemia is a frequent and serious complication affecting patient outcome. Since Microdialysis allows continuous surveillance of cerebral metabolism in a clinical setting, secondary ischemia and mitochondrial dysfunction can be recognized at an early stage. Thus, the technique opens a window of opportunity for therapeutic interventions.**

## Microdialysis sampling


Microdialysis sampling is carried out by placing the sterile FDA cleared Microdialysis catheter in the brain parenchyma. All Brain Microdialysis Catheters have a gold thread in the tip so confirmation of placement can be verified by CT.

### 70 Brain Microdialysis Catheter



- Free positioning and fixation by tunnelation
- Available in different shaft and membrane lengths

### 70 Bolt Microdialysis Catheter



- Access and fixation using a bolt system

### 106 Microdialysis Pump



The 106 Microdialysis Pump, operating at a fixed flow of 0.3 uL/min, is dedicated for the perfusion of Microdialysis catheters with sterile isotonic perfusion fluid.

## ISCUS<sup>flex</sup> Microdialysis Analyzer

The ISCUS<sup>flex</sup> Microdialysis Analyzer is specially designed to handle small sample volumes. It is a point of care analyzer for monitoring metabolic changes in brain tissue.

### Biochemical markers:

**Glucose**

**Lactate**

**Pyruvate**

**Glycerol**

**Glutamate**

**LP-ratio**



The ISCUS<sup>flex</sup> Analyzer is easily operated by medical professionals. It provides unique opportunities for early detection of metabolic crisis, ischemia and to guide post-operative interventions. Data is displayed as trend curves for easy and fast interpretation.