
## Instructions for use 70 MICRODIALYSIS BOLT CATHETER

### INTENDED USE

70 Microdialysis Bolt Catheter is intended for implantation into human brain tissue in order to enable microdialysis of the extracellular fluid of the brain.

### PRODUCT DESCRIPTION & TECHNICAL INFORMATION

X=Identical parts on all catheters

**REF  
P000131**

	<b>REF P000131</b>
1. Microvial (Polystyrene + Santoprene)	X
2. Luer lock connection (Polycarbonate)	X
3. Inlet tube (Polyurethane, OD 1 mm)	600 mm
4. Vial holder (Polycarbonate)	X
5. Outlet tube (Polyurethane, OD 1 mm)	220 mm
6. Liquid cross (ABS)	X
7. Protection tube (Polyethylene)	145 mm
8. Shaft (Polyurethane, OD 0,9 mm, Stainless steel)	130 mm
9. Dialysis membrane (Polyamide, OD 0,6 mm) Cut off: 20 000 Dalton	10 mm
10. Gold thread within the catheter membrane tip (OD 0,13 mm; L=3 mm).	X
Dead Volume Time - From membrane to microvial (at 0,3 μl/min)	23 min

The distal part of the catheter has a gold thread within the catheter tip, which makes it possible to detect the location of the catheter in the tissue through CT-scanning.

### ACCESSORIES

The 70 Microdialysis Bolt Catheter shall only be used with these accessories.

REF	Name
P000151	Perfusion fluid CNS 10x5 mL
8010191	106 Syringe 20/pkg
P000154	Microvials in rack, Sterile 12x4
P000003	106 Microdialysis Pump
P000127	107 Microdialysis Pump
8070228	MD System Complete

### CONSUMABLES

The 70 Microdialysis Bolt Catheter should be used with these consumables.

REF	Name
P000001	Microvials 250/pkg
<b>RAUMEDIC</b> , Art. 096076-001	BOLT KIT - PTO 2L, Bolt Kit for 2 brain probes
<b>INTEGRA REF IM3 Licox®</b>	IM3 Bolt Kit for 3 brain probes

## INDICATIONS

- Patients with clinical signs of brain injury or brain disease where craniotomy is required for diagnosis or therapy e.g. monitoring of ischemia in patients suffering traumatic brain injury (TBI) and subarachnoid haemorrhage (SAH). Microdialysis shall not be used as the sole basis for diagnosis or therapy.

## CONTRAINDICATIONS & RISKS

- Patients with coagulopathy, increased susceptibility to infections or bleeding disorders.
- Patients on anticoagulant drug therapy.
- When monitoring patients with brain tumors there could be a possibility of dissemination of tumor cells.
- Inserting the catheter into the brain may cause bleeding from damaged vessels.
- Leakage of cerebrospinal fluid may occur at the site of skin penetration.
- The 70 Microdialysis Bolt Catheter is for single use only. If the device is re-used there is a risk for cross-contamination.

## PRECAUTIONS

- This device is sterile unless the package has been opened or damaged.
- The 70 Microdialysis Bolt Catheter shall only be used together with the accessories described in the previous table.
- Be sure to handle the catheter carefully to avoid kinking or other damage, particularly after removal of the protection tube. Avoid contact with the dialysis membrane.
- If any visible damage is observed the catheter shall not be used.
- If there is a suspicion that the catheter has become unsterile prior to insertion the catheter shall not be used.
- The pump syringe connected to the catheter should not manually be flushed since that could damage the dialysis membrane.
- If using a 106/107 MD Pump: Check that liquid is being pumped through the catheter by inspecting the volume in the microvial **each time** the microvials are changed.  
If using MD System: Check that there are readings on the screen.
- If using 106/107 MD Pump: If there is no fluid in the collected vial, start a flush on the pump: Open the lid, wait 3 seconds and close it again. Wait for the flush (5 minutes). Check that the tubing's are not kinked and that the microvial holder needle is correctly piercing the microvial membrane. If there is still no fluid in the collected vial, the dialysis membrane might be damaged, and the catheter has to be removed.  
If using MD System: Check on screen messages.
- Remove the catheter if there is a permanent stop in the liquid flow.
- 70 Microdialysis Bolt Catheter is biocompatible up to 30 days. It may though stop working earlier because of clogging, duration of use is up to 12 days referring to literature.

## USAGE

The following procedure should be performed by a neurosurgeon under aseptic conditions.

1. Follow the instruction for inserting the bolt (see bolt manual).
2. After insertion of the bolt, verify that the compression screw is fully open.
3. Remove the protection tube by holding the liquid cross with the male luerlock connection and TURNING the protection tube counter clockwise.
4. Insert the catheter through the access for Microdialysis (INTEGRA/Temp) in the bolt (see bolt manual).
5. Fix the catheter according to the bolt instructions (see bolt manual).
6. Tighten the compression screw to stop CSF leakage.
7. Insert a microvial into the microvial holder.
8. **If you use a 106 or 107 Microdialysis Pump:** Insert a sterile microvial into the microvial holder and connect the luer lock connection to the syringe filled with CNS perfusion fluid. Remove all air bubbles in the syringe. Place the syringe in the pump and close the lid to initiate the flush. Inspect the microvial after approximately 5 minutes to see that the perfusion fluid flows through the catheter (see more detailed instructions for the pump in the pump manual).
9. **If you use the MD system:** Connect the catheter extension to the 106 Syringe filled with room temperature Perfusion Fluid CNS, make sure there are no air bubbles in the syringe. Fill up the extension by hand with Perfusion fluid and connect to the Microdialysis Catheter inlet. Place the Syringe in the MD Cartridge and follow the on-screen instructions (see more detailed instructions, in the MD System manual).

## REMOVAL OF CATHETER

The catheter is removed by gently pulling it out through the insertion site.

**NOTE:** Discontinued/removed catheters shall be handled according to the hospital routines for biohazard material.

## SYMBOL EXPLANATION & PRODUCT LIMITATIONS

	Last date of use (YYYY-MM-DD)
	LOT number
	Catalogue number
	Single use only
	See instructions for use
	Fulfils the Medical Device Directive 93/42/EEC
	Sterilised by β-radiation
	Storage temperature (4-25 °C)
	Manufacturer
	Do no use if package is damaged
Membrane	Membrane material (Polyamide, PA)
Memb. length	Membrane length (mm)
Shaft length	Shaft length (mm)
Outlet	Outlet length (mm)
Inlet	Inlet length (mm)

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