
## Instructions for use

### 71 HIGH CUT-OFF BRAIN MICRODIALYSIS CATHETER

#### **INTENDED USE**

The 71 High Cut-Off Brain Microdialysis Catheter is a single use, sterile product that is intended to monitor tissue chemistry and sample large molecules like some cytokines in the local tissue of the brain. The 71 is to be inserted by a neurosurgeon.

#### **PRODUCT DESCRIPTION & TECHNICAL INFORMATION**

X=Identical part on all catheters	REF 8010320	REF 8010331	REF 8010337
1. Microvial (Polystyrene + Santoprene)	X	X	X
2. Luer lock connection (Polycarbonate)	X	X	X
3. Inlet tube (Polyurethane, OD 1 mm)	600 mm	600 mm	600 mm
4. Microvial holder (Polycarbonate)	X	X	X
5. Outlet tube (Polyurethane, OD 1 mm)	220 mm	220 mm	220 mm
6. Stopper (Silicone)	X	X	X
7. Liquid cross (Polysulfone)	X	X	X
8. Protection tube (Polyethylene)	145 mm	145 mm	145 mm
9. Shaft (Polyurethane, OD 0,9 mm)	60 mm	60 mm	60 mm
10. Dialysis membrane(Polyarylethersulfone, OD 0,6 mm), Cut-Off: 100 000 Dalton	10 mm	20 mm	30 mm
11. Gold thread within the catheter membrane tip (OD 0,13 mm ; L=3 mm).	X	X	X
12. Tunnelating needle - accessory	X	X	X
13. Forceps - accessory	X	X	X
Dead Volume Time - From membrane to microvial (at 0,3 µl/min)	17 min	18 min	18 min

The distal part of the catheter has a gold thread within the catheter tip, which makes it possible to detect the location of the catheter in the tissue through CT-scanning.

#### **ACCESSORIES**

The 71 High Cut-Off Brain Microdialysis Catheter shall only be used with these accessories.

REF	Name
8010191	106 Syringe 20/pkg
8050151	Perfusion Fluid CNS Dextran
P000154	Microvials in rack, Sterile 12x4
P000055	Tunnelating needle
P000056	Forceps
P000003	106 Microdialysis Pump
P000127	107 Microdialysis Pump

#### **CONSUMABLES**

The 71 High Cut-Off Brain Microdialysis Catheter should be used with these consumables.

REF	Name
P000001	Microvials 250/pkg

## **INDICATIONS**

Patients with clinical signs of brain injury or brain disease where craniotomy is required for diagnosis or therapy e.g. monitoring of ischemia in patients suffering traumatic brain injury (TBI) and subarachnoid hemorrhage (SAH). Microdialysis shall not be used as the sole basis for diagnosis or therapy.

## **CONTRAINDICATIONS & RISKS**

- Patients with known hypersensitivity to Dextran.
- Patients with coagulopathy, increased susceptibility to infections or bleeding disorders.
- Patients on anticoagulant drug therapy.
- When monitoring patients with brain tumors there could be a possibility of dissemination of tumor cells.
- Inserting the catheter into the brain may cause bleeding from damaged vessels.
- Leakage of cerebrospinal fluid may occur at the site of skin penetration.
- The 71 High Cut-Off Brain Microdialysis Catheter is for single use only. If the device is re-used there is a risk for cross-contamination.

## **PRECAUTIONS**

- This device is sterile unless the package has been opened or damaged.
- The 71 High Cut-Off Brain Microdialysis Catheter shall only be used together with the accessories described in the previous table.
- Be sure to handle the catheter carefully to avoid kinking or other damage, particularly after removal of the protection tube. Avoid contact with the dialysis membrane.
- If any visible damage is observed the catheter shall not be used.
- If there is a suspicion that the catheter has become unsterile prior to insertion the catheter shall not be used.
- The pump syringe connected to the catheter should not manually be flushed since that could damage the dialysis membrane.
- Check that liquid is being pumped through the catheter by inspecting the volume in the microvial **each time** the microvials are changed.
- If there is no fluid in the collected vial, start a flush on the pump: Open the lid, wait 3 seconds and close it again. Wait for the flush (5 minutes). Check that the tubing's are not kinked and that the microvial holder needle is correctly piercing the microvial membrane. If there is still no fluid in the collected vial, the dialysis membrane might be damaged and the catheter has to be removed.
- Remove the catheter if there is a permanent stop in the liquid flow.
- 71 High Cut-Off Brain Microdialysis Catheter is biocompatible up to 30 days. It may though stop working earlier because of clogging, duration of use is up to 12 days referring to literature.

## **USAGE**

The following procedure should be performed by a neurosurgeon under aseptic conditions.

1. Tunnelate from the wound out through the scalp.
2. Insert the catheter through the tip of the tunnelating needle.
3. Withdraw the tunnelating needle.
4. Remove the protection tube by holding the liquid cross and TURNING the protection tube counter clockwise.
5. Grip the catheter SHAFT with the forceps, proximal to the membrane
6. Insert the membrane into the brain tissue through an incision in the meninges.
7. Keep the grip on the shaft while the inlet and outlet tubes are pulled straight.
8. Fix the tubing to the scalp with two sutures around the stopper.
9. Insert a sterile microvial into the microvial holder.
10. Connect the luer lock connection to the syringe filled with perfusion fluid (see instruction for the pump in use).
11. Place the syringe in the pump and close the lid to initiate the flush (see instruction as above).
12. Inspect the microvial after approximately 6 minutes to see that the perfusion fluid flows through the catheter.

## **REMOVAL OF CATHETER**

The catheter is removed by gently pulling it out through the insertion site.

**NOTE:** Discontinued/removed catheters shall be handled according to the hospital routines for biohazard material.

## **SYMBOL EXPLANATION & PRODUCT LIMITATIONS**


	Last date of use (YYYY-MM-DD)
	LOT number
	Catalogue number
	Single use only
	Manufacturer
	Consult instructions for use
	Do not use if package is damaged
	Fulfils the Medical Device Directive 93/42/EEC
	Sterilised by $\beta$ -radiation
	Storage temperature (4-25 °C)
Membrane	Membrane material (Poly(aryl)ethersulfone, P(A)ES)
Memb. length	Membrane length (mm)
Shaft length	Shaft length (mm)
Outlet	Outlet length (mm)
Inlet	Inlet length (mm)

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REF	Name
8010191	106 Syringe 20/pkg
8050151	Perfusion Fluid CNS Dextran
P000154	Microvials in Rack, Steril 12x4
P000055	Tunneling Needle
P000056	Forceps
P000003	106 Microdialysis Pump
P000127	107 Microdialysis Pump

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