MD System

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μ dialysis

Pokyny pro použití MD System

M Dialysis AB
Hammarby Fabriksväg 43
SE-120 30 Stockholm
Švédsko
E-mail: service@mdialysis.se
Tel. +46 8 470 10 20
Web: <http://www.mdialysis.se>

MD System

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
START

Perform pre-use check

- Connect the MD Amplifier cable to the MD Unit ●
- Connect MD Test Sensor to the MD Amplifier ●
- Verifying MD Amplifier ●
- Disconnect MD Test Sensor from MD Amplifier ●
- Enter the patient's ID ●

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- Insert MD Cartridge in the MD Unit ●
- Connect the MD Amplifier cable to the MD Unit ●
- Connect the MD Sensor to the MD Amplifier ●
- Place the MD Amplifier near the catheter insertion point ●
- Filling up the Catheter and Microvial. Is there liquid in the Microvial? ●
- Is there liquid in the Microvial? ●
- Connect the MD Catheter's outlet to the MD Sensor ●

NEXT

Glu
N/A

Lac
N/A

Pyr
N/A

L/P
N/A

HOME
VIEW
SETTINGS
SERVICE
PATIENT DATA

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

START

Perform pre-use check

Connect the MD Amplifier cable to the MD Unit ●

Connect MD Test Sensor to the MD Amplifier ●

Verifying MD Amplifier ●

Disconnect MD Test Sensor from MD Amplifier ●

Enter the patient's ID ●

Fill Syringe Orange with 2.5 mL of MD Calibrator and connect to the MD Sensor

Let the Perfusion Fluid reach room temperature. Fill the 106 Syringe (non-colored) with approx 2.7 mL of sterile Perfusion Fluid.

Connect the 106 Syringe to the Catheter Extension and fill up the Catheter Extension by hand until there's 2.5 mL Perfusion Fluid left in the 106 Syringe

Connect the filled up Catheter Extension with the inlet of the Microdialysis catheter. Put a Microvial in the Microvial Holder of the Catheter Outlet.

Insert syringes into MD Cartridge. Syringe Orange in pos. 1 (Left), 106 Syringe in pos. 3 (Right)

Insert MD Cartridge in the MD Unit

Connect the MD Amplifier cable to the MD Unit ●

Connect the MD Sensor to the MD Amplifier ●

Place the MD Amplifier near the catheter insertion point ●

Filling up the Catheter and Microvial, is there liquid in the Microvial? ●

Is there liquid in the Microvial? ●

Connect the MD Catheter's outlet to the MD Sensor

Glu
N/A

Lac
N/A

Pyr
N/A

L/P
N/A

HOME
VIEW
SETTINGS
SERVICE
PATIENT DATA
START

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START

Perform pre-use check

Connect the MD Amplifier cable to the MD Unit ●

Connect MD Test Sensor to the MD Amplifier ●

Verifying MD Amplifier ●

Disconnect MD Test Sensor from MD Amplifier ●

Enter the patient's ID ●

Fill Syringe Orange with 2.5 mL of MD Calibrator and connect to the MD Sensor

Let the Perfusion Fluid reach room temperature. Fill the 106 Syringe (non-colored) with approx 2.7 mL of sterile Perfusion Fluid.

Connect the 106 Syringe to the Catheter Extension and fill up the Catheter Extension by hand until there's 2.5 mL Perfusion Fluid left in the 106 Syringe

Connect the filled up Catheter Extension with the inlet of the Microdialysis catheter. Put a Microvial in the Microvial Holder of the Catheter Outlet.

Insert syringes into MD Cartridge. Syringe Orange in pos. 1 (Left), 106 Syringe in pos. 3 (Right)

Insert MD Cartridge in the MD Unit ●

Connect the MD Amplifier cable to the MD Unit ●

Connect the MD Sensor to the MD Amplifier ●

Place the MD Amplifier near the catheter insertion point ●

Filling up the Catheter and Microvial, is there liquid in the Microvial?

Is there liquid in the Microvial? ●

Connect the MD Catheter's outlet to the MD Sensor

YES

Glu
N/A

Lac
N/A

Pyr
N/A

L/P
N/A

HOME
VIEW
SETTINGS
SERVICE
PATIENT DATA
START

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The screenshot shows the MD Monitor interface during a 'PAUSE' state. The screen displays a diagram of the MD Cartridge being ejected from the MD Unit, with a 'CLICK' label and an arrow indicating the direction of movement. The background shows the same real-time data as the previous screenshot: Glu (4,8), Lac (2,2), Pyr (90,7), and L/P (24,4). The interface includes a 'RESUME' button and various status indicators like 'MICROVIAL SAMPLE', '2021-04-22', and '12:31:42'.

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
PAUSE

Retracting motors and ejecting cartridge

Disconnect the catheter vial connector from the MD Sensor

Insert Microvial in Microvial holder

Remove MD Cartridge

Remove the 106 syringe from the MD Cartridge

Remove the MD Amplifier with the MD Sensor from the patient

Place the MD Amplifier near the Catheter insertion point

Remove the Microvial and connect the MD Catheter's outlet to the MD Sensor

Reinsert 106 syringe in the MD Cartridge

Reload the MD Cartridge in the MD Unit

Glu
N/A

Lac
N/A


Pyr
N/A

L/P
N/A

NEXT

VIEW SETTINGS SERVICE PATIENT DATA RESUME MICROVIAL SAMPLE 2021-04-22 11:30:45 LOKE μdialysis

6. Odpojte MD Amplifier s MD Sensor em od pacienta.

PAUSE

Retracting motors and ejecting cartridge

Disconnect the catheter vial connector from the MD Sensor

Insert Microvial in Microvial holder

Remove MD Cartridge

Remove the 106 syringe from the MD Cartridge

Remove the MD Amplifier with the MD Sensor from the patient

Place the MD Amplifier near the Catheter insertion point

Remove the Microvial and connect the MD Catheter's outlet to the MD Sensor

Reinsert 106 syringe in the MD Cartridge

Reload the MD Cartridge in the MD Unit

Glu
N/A

Lac
N/A


Pyr
N/A

L/P
N/A

NEXT

VIEW SETTINGS SERVICE PATIENT DATA RESUME MICROVIAL SAMPLE 2021-04-22 11:33:32 LOKE μdialysis

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The screenshot displays the control panel in a 'PAUSE' state. The main instruction area lists the following steps:

- Retracting motors and ejecting cartridge
- Disconnect the catheter vial connector from the MD Sensor
- Insert Microvial in Microvial holder
- Remove MD Cartridge
- Remove the 106 syringe from the MD Cartridge
- Remove the MD Amplifier with the MD Sensor from the patient
- Place the MD Amplifier near the Catheter insertion point**
- Remove the Microvial and connect the MD Catheter's outlet to the MD Sensor
- Reinsert 106 syringe in the MD Cartridge
- Reload the MD Cartridge in the MD Unit

On the right side, the following values are displayed:

- Glu: N/A
- Lac: N/A
- Pyr: N/A
- L/P: N/A

The bottom navigation bar includes icons for Home, View, Settings, Service, Patient Data, a highlighted Resume button, Microvial Sample, the date and time (2021-04-22 11:35:24), a battery icon, and the LOKE μdialysis logo.

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The screenshot displays the control panel in a 'PAUSE' state. The main instruction area lists the following steps:

- Retracting motors and ejecting cartridge
- Disconnect the catheter vial connector from the MD Sensor
- Insert Microvial in Microvial holder
- Remove MD Cartridge
- Remove the 106 syringe from the MD Cartridge
- Remove the MD Amplifier with the MD Sensor from the patient
- Place the MD Amplifier near the Catheter insertion point
- Remove the Microvial and connect the MD Catheter's outlet to the MD Sensor**
- Reinsert 106 syringe in the MD Cartridge
- Reload the MD Cartridge in the MD Unit


On the right side, the following values are displayed:

- Glu: N/A
- Lac: N/A
- Pyr: N/A
- L/P: N/A

The bottom navigation bar includes icons for Home, View, Settings, Service, Patient Data, a highlighted Resume button, Microvial Sample, the date and time (2021-04-22 11:36:01), a battery icon, and the LOKE μdialysis logo.

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PAUSE

Retracting motors and ejecting cartridge

Disconnect the catheter vial connector from the MD Sensor
Insert Microvial in Microvial holder
Remove MD Cartridge
Remove the 106 syringe from the MD Cartridge
Remove the MD Amplifier with the MD Sensor from the patient
Place the MD Amplifier near the Catheter insertion point
Remove the Microvial and connect the MD Catheter's outlet to the MD Sensor

Reinsert 106 syringe in the MD Cartridge

Reload the MD Cartridge in the MD Unit

Glu
N/A

Lac
N/A

Pyr
N/A


L/P
N/A

NEXT

HOME VIEW SETTINGS SERVICE PATIENT DATA RESUME MICROVIAL SAMPLE 2021-04-22
11:36:47 LOKE
μdialysis

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PAUSE

Retracting motors and ejecting cartridge

Disconnect the catheter vial connector from the MD Sensor
Insert Microvial in Microvial holder
Remove MD Cartridge
Remove the 106 syringe from the MD Cartridge
Remove the MD Amplifier with the MD Sensor from the patient
Place the MD Amplifier near the Catheter insertion point
Remove the Microvial and connect the MD Catheter's outlet to the MD Sensor
Reinsert 106 syringe in the MD Cartridge

Reload the MD Cartridge in the MD Unit

Glu
N/A

Lac
N/A

Pyr
N/A

L/P
N/A

RUN

HOME VIEW SETTINGS SERVICE PATIENT DATA RESUME MICROVIAL SAMPLE 2021-04-22
11:41:36 LOKE
μdialysis

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

5.1. Error list

#	Text
-1	-1: "Unknown error",
0	0: "ADC Process: Unable to initialize ADC after {retries} attempts, disabling ADC.",
1	1: "Unable to connect to database - {e}.",
2	2: "Driver/Controller initialize error - {e}.",
3	3: "Waiting for processes to load timed out, exiting.",
4	4: "Cartridge Controller Eject: Error ejecting cartridge - {e}.",
5	5: "LED Controller LEDS1 Off: Error turning off LEDs - {e}.",
6	6: "LED Controller LEDS2 Off: Error turning off LEDs - {e}.",
7	7: "{module_name}: Calibration Error: Calibration factor needs to be cleared before gathering measurements for new calculation.",
8	8: "{module_name}: Calibration Error: No measurements stored.",
9	9: "{module_name}: Average signal from biosensor is 0, cannot calculate calibration factor.",
10	10: "{module_name}: Average signal from biosensor is less than 0 ({avg_snorm}), cannot calculate calibration factor.",
11	11: "{module_name}: Factory calibration factor is 0, cannot do comparison.",
12	12: "{module_name}: New calibration factor {calibration_factor} differs more than {cal_factory_cal_max_difference} from factory calibration: {factory_cal}.",
13	13: "Unable to parse test script file: {e}.",
14	14: "Patient Module Controller Biosensor Detect Callback: Could not read biosensor EEPROM, type of biosensor could not be determined.",
15	15: "Patient Module Controller Read EEPROM Data: {e}.",
16	16: "Patient Module Controller Measure Temperature - Could not measure temperature after {retries} retries: {e}.",
17	17: "Perfusion Flow rate not defined in Analysis Controller, cannot perform flush - aborting analysis.",
18	18: "A pump task for perfusion has already been defined in the script, cannot perform flush - aborting analysis.",
19	19: "Flush: Error in motor task list - aborting analysis.",
20	20: "Biosensor removed during flush, aborting analysis.",
21	21: "Cartridge forcibly ejected during flush, aborting analysis.",
22	22: "Unable to write start time to EEPROM: {e}.",
23	23: "Perfusion syringe empty. Re-fill syringe.",
24	24: "Error during pumping, aborting analysis.",
25	25: "Analysis Error - Biosensor expired on {expired_dt_str}.",
26	26: "Biosensor removed, aborting analysis.",
27	27: "Cartridge forcibly ejected, aborting analysis.",
28	28: "Sensor temperature sensor error, replace sensor.",
29	29: "Sensor Error - Short circuit detected, replace sensor.",
30	30: "Calibrator syringe empty. Re-fill syringe.",
31	31: "Calibration failed after {attempts} attempts. Check syringe placement or replace sensor.",
32	32: "Biosensor removed during calibration, aborting analysis.",
33	33: "Cartridge forcibly ejected during calibration, aborting analysis.",
34	34: "Invalid patient id/name ({patient_id} - {patient_name}) or flow rate {flow_rate}, aborting analysis.",
35	35: "Analysis Error - Error in motor task list during fill-up script.",
36	36: "Analysis Error - Timeout during fill up.",
37	37: "System error. Check syringes or replace sensor.",
38	38: "Analysis Error - Error in motor task list during fill up post script.",

#	Text
39	39: "Biosensor removed during fill-up sequence, aborting analysis.",
40	40: "Cartridge forcibly ejected during fill-up sequence, aborting analysis.",
41	41: "Biosensor removed, aborting start of analysis.",
42	42: "Cartridge forcibly ejected, aborting analysis.",
43	43: "Find Piston Error - Error in motor task list during find piston.",
44	44: "Timed out during Pre-Use Check while waiting for measurement data.",
45	45: "Biosensor Channel {channel} limit test failed, value: {ncurrent}.",
46	46: "Could not fully retract motors, remove cartridge carefully by hand and contact service.",
47	47: "Date/Time for MD Unit set incorrectly ({timestamp}), contact service.",
48	48: "Error during initialization of motors, contact service.",
49	49: "Flush Error - Error in motor task list.",
50	50: "Biosensor removed during flush, aborting analysis.",
51	51: "Cartridge forcibly ejected during flush, aborting analysis.",
52	52: "Could not fully retract motors, remove cartridge carefully by hand.",
53	53: "Could not fully retract motors.",
54	54: "Could not fully retract motors, aborting test mode.",
55	55: "Error in task list ID {task_list_id}: {task_list_errors}.",
56	56: "Error in task {task_motor_name} / {task_task_name}: {task_errors}.",
57	57: "Error in battery info thread: {e}.",
58	58: "Biosensor data signal missing check failed.",
59	59: "{key}: Signal missing - max value last {duration}s: {ncurrent} nA below limit {limit_na} nA.",
60	60: "{key}: Signal too high - average value last {duration}s: {ncurrent} nA over limit {limit_na} nA.",
61	61: "Biosensor data short circuit check failed.",
62	62: "Temperature below allowed low limit ({limit_low}) - average value last {temperature_sanity_length} measurements: {avg_temp:.2f}°C.",
63	63: "Temperature above allowed high limit ({limit_high}) - average value last {temperature_sanity_length} measurements: {avg_temp:.2f}°C.",
64	64: "Unable to get channel average ncurrent: {e}.",
65	65: "Resetting MD Amplifier by power cycle...",
66	66: "Already tried resetting MD Amplifier once, temperature cannot be read.",
67	67: "Pre-Use Check Failed, values: {channel_data}.",
68	68: "Error quering the max date: {e}.",
69	69: "Error inserting into database: {e}.",
70	70: "Unable to parse request data: {json_data}, error: {e}.",
71	71: "Unknown LogMessage level: {log_message}.",
72	72: "Unable to parse push data: {data}, error: {e}.",
73	73: "Socket Event Handler Exception: {e}, Event Data: {event_data}",
74	74: "Invalid socket request received: {request}.",
75	75: "Socket Request Handler Exception: {e}.",
76	76: "ADS1254 Driver: Unable to read at least 1 sample.",
77	77: "ADS1254 Driver: Max read sample attempts reached.",
78	78: "ADS1254 Driver: Error reading from ADC - {e}.",
79	79: "ADS1254 Driver: Read bits timeout, couldn't find DREADY signal pattern.",
80	80: "Exception: {e}.",
81	81: "Unable to merge user settings with default settings: {e}.",
82	82: "Could not copy fallback file to destination: {e}.",
83	83: "Unable to merge fallback user settings with default settings: {e}.",

#	Text
84	84: "Unable to merge user settings debug file with default settings: {e}."
85	85: "Unable to load settings file: {e}."
86	86: "Could not find path {element} in settings."
87	87: "ADC Process: Error intializing ADC, retry count {init_retry_count} - {e}."
88	88: "ADC Process: Unable to set ADC Registry Configuration, retry count {init_retry_count}."
89	89: "Motor Process: Unable to initialize Serial Port Driver - {e}."
90	90: "Possible air bubble detected - Difference between measurement type {mtype} and {reftype} greater than {max_diff}, minimum difference last {min_duration}s = {min_diff}. Current {reftype} measurement value: {ref_meas}mmol/L."
91	91: "Unable to calculate biosensor reference levels - {e}."
92	92: "{key}: Signal differs too much from reference level - max value last {duration}s: {ncurrent} nA over limit {limit_na} nA."
93	93: "Reference level check has failed."
94	94: "Sensor Error - Signal missing error detected, replace sensor."
95	95: "Sensor Error - Possible air bubble detected, unable to automatically recover."
96	96: "Sensor Error - Reference level error possibly due to stop in flow, unable to automatically recover."
97	97: "{key}: Signal too high - min value last {duration}s: {ncurrent} nA over limit {limit_na} nA."
98	98: "Post calibration check has failed."
99	99: "Sensor Error - Error possibly due to stop in flow, unable to automatically recover."
100	100: "Post calibration check failed for {name} - min value last {duration}s: {min_value} mmol/L over limit {limit} mmol/L."
101	101: "Aborting analysis after {attempts} attempts at filling up the Microvial."
102	102: "Take Sample duration invalid: {sample_duration}."

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