REAGENT Set C

For ISCUS/ISCUS^{flex} Clinical Microdialysis Analyzer

Ref. No. 8002165

Intended user: Medical or laboratory professional staff
Content:

1. Reagent: One bottle lyophilized reagent each for Glucose, Lactate and Pyruvate placed in a cassette.
2. Buffer: One bottle 6 mL each for Glucose, Lactate and Pyruvate
3. Calibrator: One bottle of Calibrator A 6 mL placed in the Reagent Cassette.

Preparation

1. Unscrew the cap with the membrane from the reagent and calibrator bottles, placed in the cassette. Remove and discard the rubber stoppers.
2. Unscrew the cap from the buffer bottles and gently transfer the contents to the corresponding reagent bottle.
3. Fasten the cap with the membrane on the reagent and calibrator bottles in the cassette, without returning the rubber stoppers.
4. Let the reagents and calibrator stand and equilibrate in room temperature for at least 30 minutes prior to use. The contents will be completely mixed when the Reagent Cassette is placed and identified in the analyzer.

Intended purpose and stability of solution


he Reagent Set C is a cassette containing reagents and Calibrator A made for use in ISCUS/ISCUS^{flex} Microdialysis Analyzer. For intended purpose, see information for the individual components. Each Reagent Cassette has a unique code, written on the label, which should be registered when placing it in the analyzer. The Reagent Set is stable up to expiry date when stored at +2 to +8 °C. Reconstituted reagents are stable for five days in the analyzer. The contents are sufficient for around 310 analyses.

GLUCOSE

Colorimetric method for the quantitative determination of Glucose in Microdialysates.

Measuring principle

Glucose is enzymatically oxidized by glucose oxidase (GOD). The hydrogen peroxide formed reacts with phenol and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields a red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the glucose concentration.

Linear range: 0.1 - 25 mmol/L

	Component	Concentration in test solution
Glucose reagent	4-aminoantipyrine	0.77 mmol/L
	Ascorbate oxidase	>3 kU/L
	Glucose oxidase	>1.5 kU/L
	Peroxidase	>1.5 kU/L
Glucose buffer	Phosphate buffer, pH 7.0	0.1 mmol/L
	Phenol	11 mmol/L
	Sodium azide	0.4 g/L

References:



1. Barhem and P. Trinder, Analyst 97 (1972) 142

LACTATE

Colorimetric method for the quantitative determination of Lactate in Microdialysates.

Measuring principle

Lactate is enzymatically oxidized by lactate oxidase. The hydrogen peroxide formed reacts with 4-chlorophenol and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields a red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the lactate concentration.

Linear range: 0.1 - 12 mmol/L

	Component	Concentration in test solution
Lactate reagent	4-aminoantipyrine	0.4 mmol/L
	Lactate oxidase	>500 U/L
	Peroxidase	>500 U/L
	Ascorbate oxidase	>12.0 kU/L
Lactate buffer	PIPES buffer, pH 6.8	100 mmol/L
	4-Chlorophenol	5.4 mmol/L
	Sodium oxalate	7.5 mmol/L
	EDTA-disodium salt	5 mmol/L
	Sodium azide	0.3 g/L

References:


1. N. Shimojo et al. Clin Chem 35 (1989) 1992
2. H.F. Kühnle et al. J.Clin Chem BioChem 15 (1977) 171
3. T.O. Kleine et al. Dtsch Med Wschr 104 (1979) 553

PYRUVATE

Colorimetric method for the quantitative determination of Pyruvate in Microdialysates.

Measuring principle

Pyruvate is enzymatically oxidized by pyruvate oxidase (PyrOx). The hydrogen peroxide formed reacts with N-ethyl-N-(2-hydroxy-3-sulfo-propyl)-m-toluidine and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields a red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the pyruvate concentration.

Default Linear range: 10 - 300 µmol/L

	Component	Concentration in test solution
Pyruvate reagent	4-amino-antipyrine	0.3 mmol/L
	Tiamine pyrophosphate	0.2 mmol/L
	FAD	10 µmol/L
	Pyruvate Oxidase	>0.2 kU/L
	Peroxidase	>0.8 kU/L
	Ascorbate Oxidase	>10 kU/L
Pyruvate buffer	Citrate buffer, pH 6.1	100 mmol/L
	Potassium dihydrogenphosphate	10 mmol/L
	MgCl2	10 mmol/L
	TOOS	1.5 mmol/L
	Sodium azide	0.3 g/L

References:

1. B. Sedewitz, et al., J. Bacteriol, 160 (1984) 273-278
2. M. Nawata, et al., Anal Biochem, 190 (1990) 84-87
3. H. Araki and M. Yamada, in: H. U. Bergmeyer (Editor), Methods of Enzymatic Analysis, 3rd ed., Vol 6, Verlag Chemie, Weinheim, 1984

CALIBRATOR A

Calibration values	
Glucose	5.55 mmol/L
Lactate	2.5 mmol/L
Pyruvate	250 µmol/L

Sample material

Microdialysates

For in vitro use only

Symbol declaration

Last day of use

Lot number

Storage temperature

Consult instructions for use

In vitro diagnostic device

The product meets EU directive for IVD (98/79/EC) LVFS 2001:7

Manufactured by:
M Dialysis AB
Hammarby Fabriksväg 43
SE-120 30 · Stockholm · Sweden
Tel: +46-8-470 10 20
Fax: +46-8-470 10 55
E-mail: info@mdialysis.com
www.mdialysis.com

USA office:
M Dialysis Inc.
73 Princeton Street
N.Chelmsford · MA 01863 · USA
Phone: +1 978 251 1940, +1-866-868-9236
Fax: +1 978 251-1960
E-mail: usa@mdialysis.com

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