Reagent Kit


for the 600 & ISCUS^{flex}

Microdialysis Analyzer

GLUCOSE

Intended purpose: Colorimetric method for the quantitative determination of Glucose in Microdialysates.

Measuring principle
Glucose is enzymatically oxidised by glucose oxidase (GOD). The hydrogen peroxide formed reacts with phenol and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields the red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the glucose concentration.

Linear range: 0.1 - 25 mmol/L

	Component	Concentration in test solution
Glucose reagent	4-Aminoantipyrine	0.77 mmol/L
	Ascorbate oxidase	>3 kU/L
	Glucose oxidase	>1.5 kU/L
	Peroxidase	>1.5 kU/L
Glucose buffer	Phosphate buffer, pH 7.0	0.1 mol/L
	Phenol	11 mmol/L
	Sodium azide	0.4 g/L

Sample material

Microdialysates

WARNING:

Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.
The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10 % sodium hydroxide

Symbol declaration:

Last day of use

LOT

Lot number

Storage temperature

Consult instructions for use

IVD

In vitro diagnostic device

CE

The product meets EU directive for IVD (98/79/EC) LVFS 2001:7

For in vitro use only

References:
1. Barhem and P.Trinder, Analyst 97(1972)142

Reagent Kit


for the 600 & ISCUS^{flex}

Microdialysis Analyzer

LACTATE

Intended purpose: Colorimetric method for the quantitative determination of Lactate in Microdialysates.

Measuring principle
Lactate is enzymatically oxidised by lactate oxidase. The hydrogen peroxide formed reacts with 4-chlorophenol and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields the red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the lactate concentration.

Linear range: 0.1 - 12 mmol/L

	Component	Concentration in test solution
Lactate reagent	4-Aminoantipyrine	0.4 mmol/L
	Lactate oxidase	>0.5 kU/L
	Peroxidase	>0.5 kU/L
	Ascorbate oxidase	>12.0 kU/L
Lactate buffer	PIPES buffer, pH 6.8	100 mmol/L
	4-Chlorophenol	5.4 mmol/L
	Sodium oxalate	7.5 mmol/L
	EDTA-disodium salt	5 mmol/L
	Sodium azide	0.3 g/L

Sample material

Microdialysates

WARNING:

Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.
The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10 % sodium hydroxide

Symbol declaration:

Last day of use

LOT

Lot number

Storage temperature

Consult instructions for use

IVD

In vitro diagnostic device

CE

The product meets EU directive for IVD (98/79/EC) LVFS 2001:7

For in vitro use only

References:
1.N.Shimojo et al. Clin Chem 35(1989)1992 2.H.F.Kühnle et al, J.Clin Chem BioChem 15 (1977)171
2.T.O.Kleine et al, Dtsch Med Wschr 104 (1979) 553

Reagent Kit


for the 600 & ISCUS^{flex}

Microdialysis Analyzer

PYRUVATE

Intended purpose: Colorimetric method for the quantitative determination of Pyruvate in Microdialysates.

Measuring principle
Pyruvate is enzymatically oxidized by pyruvate oxidase (PyrOx). The hydrogen peroxide formed reacts with N-ethyl-N-(2-hydroxy-3-sulfopropyl)-m-toluidine and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields the red-violet colored quinonediiimine. The rate of formation is measured photometrically at 530 nm and is proportional to the pyruvate concentration.

Deafault Linear range: 10 - 300 µmol/L

	Component	Concentration in test solution
Pyruvate reagent	4-Aminoantipyrine	0.3 mmol/L
	Tiamine pyrophosphate	0.2 mmol/L
	FAD	10 µmol/L
	Pyruvate Oxidase	>0.2 kU/L
	Peroxidase	>0.8 kU/L
Pyruvate buffer	Ascorbate Oxidase	>10 kU/L
	Citrate buffer, pH 6.1	100 mmol/L
	Potassium dihydrogenphosphate	10 mmol/L
	MgCl ₂	10 mmol/L
	TOOS	1.5 mmol/L
	Sodium azide	0.3 g/L

Sample material

Microdialysates

WARNING:

Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.
The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10 % sodium hydroxide

Symbol declaration:

Last day of use

LOT

Lot number

Storage temperature

Consult instructions for use

IVD

In vitro diagnostic device

CE

The product meets EU directive for IVD (98/79/EC) LVFS 2001:7

For in vitro use only

References:
1. B. Sedewitz, et al., J. Bacteriol., 160 (1984) 273-278
2. M. Nawata, et al., Anal Biochem., 190 (1990) 84-87
3. H. Araki and M. Yamada, in: H. U. Bergmeyer (Editor), Methods of Enzymatic Analysis, 3rd ed., Vol 6, Verlag Chemie, Weinheim, 1984

Reagent Kit

for the 600 & ISCU^{flex}


Microdialysis Analyzer

GLYCEROL

Intended purpose: Colorimetric method for the quantitative determination of Glycerol in Microdialysates.

Measuring principle

Glycerol is phosphorylated by adenosine triphosphate (ATP) and glycerol kinase (GK) to glycerol-3-phosphate, which is subsequently oxidized in the presence of glycerol-3-phosphate oxidase (GPO). The hydrogen peroxide formed reacts with 3,5-dichloro-2-hydroxy-benzene sulphonic acid (DCHBS) and 4-amino-antipyrine. This reaction is catalyzed by peroxidase

Linear range: 0.01 - 1.5 mmol/L

	Component	Concentration in test solution
Glycerol reagent	4-aminoantipyrine	0.4 mmol/L
	ATP	1.0 mmol/L
	Glycerol kinase	>400 U/L
	Glycerol-3-phosphate-oxidase	>1.5 kU/L
	Peroxidase	>1 kU/L
Glycerol buffer	Ascorbate oxidase	>7.0 kU/L
	PIPES buffer, pH 7.6	40 mmol/L
	DCHBS	1.5 mmol/L
	Magnesium ions	17.5 mmol/L

Sample material
Microdialysates

Symbol declaration:

 Last day of use

 Lot number

 Storage temperature

 Consult instructions for use

 In vitro diagnostic device

 The product meets EU directive for IVD (98/79/EC) LVFS 2001:7

WARNING:

Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.

The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.

Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10 % sodium hydroxide

For in vitro use only

Content


REF P000011

GB

1. Reagent: One bottle lyophilized reagent each for glucose, lactate, pyruvate and glycerol.
 2. Buffer: One bottle à 6 mL each for glucose, lactate, pyruvate and glycerol.
 3. Calibrator: One bottle à 6 mL.
- Reagents are sufficient for 310 determinations.
Reagents and calibrator are stable up to expiry date when stored at +2 to +8 °C.

Preparation and stability of solution

1. Unscrew the cap with the membrane from the reagent bottle. Remove and discard the rubber stopper.
 2. Transfer the contents of the buffer bottle to the reagent bottle.
 3. Fasten the cap with the membrane on the reagent bottle, without Rubber stopper.
 4. Dissolve contents completely by gently turning the bottle upside-down at least ten times. Let the reagent stand and equilibrate in room temperature for at least 30 minutes prior to use.
- Reconstituted reagent is stable for five days in the instrument.

- Dissolve contents completely by gently turning the bottle upside-down at least ten times.

Intended user: Medical or laboratory professional staff

Intended purpose: see information for the individual components.

Manufactured by:
M Dialysis AB
Hammarby Fabriksväg 43
SE-120 30 • Stockholm • Sweden
Tel: +46-8-470 10 20
Fax: +46-8-470 10 55
E-mail: info@mdialysis.com
www.mdialysis.com

USA office:
M Dialysis Inc.
73 Princeton Street
N. Chelmsford • MA 01863 • USA
Tel: (800) 440-4980, (978) 251-1940
Fax: (978) 251-1960
E-mail: usa@mdialysis.com

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
Referenser:
1. Barhem and P.Trinder, Analyst 97(1972)142

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