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

for the 600 and ISCUS^{flex}
Microdialysis Analyzers

LACTATE

REF. P000024
Intended purpose: Colorimetric method for the quantitative determination of lactate in Microdialysates.
Also required: Calibrator A. REF. P000057 & Rinsing Fluid, REF. 8002171

Measuring principle
Lactate is enzymatically oxidised by lactate oxidase (LOD). The hydrogen peroxide formed reacts with 4-chlorophenol and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields the red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the lactate concentration.

- Reagents:
1. Reagent: 5 bottles of lyophilisate
 2. Buffer: 5 bottles of 6 mL
- Reagent sufficient for 5 x 350 determinations.
Reagents are stable up to expiry date when stored at +2 to +8°C.

- Preparation and stability of solution
1. Unscrew the cap with the membrane from the reagent bottle. Remove and discard the rubber stopper.
 2. Transfer the contents of the buffer bottle to the reagent bottle.
 3. Fasten the cap with the membrane on the reagent bottle, without Rubber stopper.
 4. Dissolve contents completely by gently turning the bottle upside-down at least ten times. Let the reagent stand and equilibrate in room temperature for at least 30 minutes prior to use.
- Reconstituted reagent is stable for five days in the instrument.

	Component	Concentration in test solution
Lactate reagent	4-aminoantipyrine	0.4 mmol/L
	Lactate oxidase	>500 U/L
	Peroxidase	>500 U/L
Lactate buffer	Ascorbate oxidase	>12.0 kU/L
	PIPES buffer, pH 6.8	100 mmol/L
	4-Chlorophenol	5.4 mmol/L
	Sodium oxalate	7.5 mmol/L
	EDTA-disodium salt	5 mmol/L
	Sodium azide	0.3 g/L

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See instructions for use

In vitro diagnostic reagent

The product meets EU directive for IVD (98/79/EC)/LVFS 2001:7

Calibration

Use Calibrator A. REF. P000057

WARNING:

Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.

The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.

Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10 % sodium hydroxide

Referenser: 1. N.Shimojo et al. Clin Chem 35(1989)1992 2.H.F.Kühnle et al, J.Clin Chem BioChem 15 (1977)171. 3. T.O.Kleine et al, Dtsch Med Wscht 104 (1979) 553.

References: 1. N.Shimojo et al. Clin Chem 35(1989)1992 2.H.F.Kühnle et al, J.Clin Chem BioChem 15 (1977)171. 3. T.O.Kleine et al, Dtsch Med Wscht 104 (1979) 553.

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GB: Quality Control:

It is recommended that the system be controlled using Control Samples. These controls should be assayed according to local quality assurance schemes. The use of other control materials has not been evaluated. Satisfactory level of performance is achieved when the analyte values for the controls are within the "Acceptable Control Range" published in the Packet Insert with the controls.

Additional information is found in the Technical Manual for the Microdialysis Analyzers.

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M Dialysis AB