REAGENS

SE

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GLYCEROL

REF: P000025

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GLYCEROL


REF: P000025

Intended purpose: Colorimetric method for the quantitative determination of glycerol in Microdialysates.

Also required: Calibrator A. REF: P000057 & Rinsing Fluid, REF. 8002171

Measuring principle

Glycerol is phosphorylated by adenosine triphosphate (ATP) and glycerol kinase (GK) to glycerol-3-phosphate, which is subsequently oxidized in the presence of glycerol-3-phosphate oxidase (GPO). The hydrogen peroxide formed reacts with 3,5-dichloro-2-hydroxy-benzene sulphonic acid (DCHBS) and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields the red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the glycerol concentration.

Linear range: 0.01 - 1.5 mmol/L

Reagents:

1. Reagent: 5 bottles of lyophilisate
 2. Buffer: 5 bottles of 6 mL
- Reagent sufficient for 5 x 350 determinations.
Reagents are stable up to expiry date when stored at +2 to +8°C.

Preparation and stability of solution

1. Unscrew the cap with the membrane from the reagent bottle. Remove and discard the rubber stopper.
 2. Transfer the contents of the buffer bottle to the reagent bottle.
 3. Fasten the cap with the membrane on the reagent bottle, without Rubber stopper.
 4. Dissolve contents completely by gently turning the bottle upside-down at least ten times. Let the reagent stand and equilibrate in room temperature for at least 30 minutes prior to use.
- Reconstituted reagent is stable for five days in the instrument.

	Component	Concentration in test solution
Glycerol reagent	4-aminoantipyrine	0.4 mmol/L
	ATP	1.0 mmol/L
	Glycerol kinase	>400 U/L
	Glycerol-3-phosphate-oxidase	>1.5 kU/L
	Peroxidase	>1 kU/L
Glycerol buffer	Ascorbate oxidase	>7.0 kU/L
	PIPES buffer, pH 7.6	40 mmol/L
	DCHBS	1.5 mmol/L
	Magnesium ions	17.5 mmol/L
	Sodium azide	0.2 g/L

Sample material
Microdialysates

For in vitro use only

Symbol declaration:

-  Last day of use
-  Lot number
-  Storage temperature
-  See instructions for use
-  In vitro diagnostic reagent
-  The product meets EU directive for IVD (98/79/EC)/LVFS 2001:7

Calibration
Use Calibrator A. REF. P000057

WARNING:

Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.

The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.

Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide

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GB: Quality Control:

It is recommended that the system be controlled using Control Samples. These controls should be assayed according to local quality assurance schemes. The use of other control materials has not been evaluated. Satisfactory level of performance is achieved when the analyte values for the controls are within the "Acceptable Control Range" published in the Packet Insert with the controls.

Additional information is found in the Technical Manual for the Microdialysis Analyzers.

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M Dialysis AB