REAGENS

SE

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
UREA

REF. P000026

Intended purpose: UV-method for the quantitative determination of urea in Microdialysates.
Also required: Calibrator A. REF: P000057 & Rinsing Fluid, REF. 8002171

Measuring principle

Urea is hydrolyzed in the presence of urease to ammonium ions and carbon dioxide. The ammonium ions react with 2-oxoglutarate in the presence of glutamate dehydrogenase (GIDH) and NADH to form glutamate and NAD⁺. The rate of utilization of NADH is measured photometrically at 375 nm and is proportional to the urea concentration.

Linear range: 0.5 - 25 mmol/L (recently prepared), 0.5 - 17 mmol/L (after three days)

Reagents

1. Reagent: 5 bottles of lyophilisate
2. Buffer: 5 bottles of 6 mL.

Reagent sufficient for 5 x 350 determinations.

Reagents are stable up to expiry date when stored at +2 to +8°C.

Preparation and stability of solution

1. Unscrew the cap with the membrane from the reagent bottle. Remove and discard the rubber stopper.
 2. Transfer the contents of the buffer bottle to the reagent bottle.
 3. Fasten the cap with the membrane on the reagent bottle, without Rubber stopper.
 4. Dissolve contents completely by gently turning the bottle upside-down at least ten times. Let the reagent stand and equilibrate in room temperature for at least 30 minutes prior to use.
- Reconstituted reagent is stable for three days in the instrument.

	Component	Concentration in test solution
Urea reagent	Urease	>15 kU/L
	Glutamate dehydrogenase	>1 kU/L
	NADH	0.2 mmol/L
	2-oxoglutarate	15 mmol/L
	TRIS buffer, pH 8.0	150 mmol/L
Urea buffer	Sodium azide	0.3 g/L

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Sample material Microdialysates	Calibration Use Calibrator A. REF. P000057
For in vitro use only	WARNING: Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents. The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention. Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide.
Symbol declaration:	
Last day of use	
Lot number	
Storage temperature	
See instructions for use	
In vitro diagnostic reagent	
The product meets EU directive for IVD (98/79/EC)/LVFS 2001:7	

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
GB: Quality Control:

It is recommended that the system be controlled using Control Samples. These controls should be assayed according to local quality assurance schemes. The use of other control materials has not been evaluated. Satisfactory level of performance is achieved when the analyte values for the controls are within the "Acceptable Control Range" published in the Packet Insert with the controls.

Additional information is found in the Technical Manual for the Microdialysis Analyzers.

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M Dialysis AB