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## Microdialysis Analyzers

## PYRUVATE

REF. P000063

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## PYRUVATE


REF. P000063

Intended purpose: Colorimetric method for the quantitative determination of pyruvate in Microdialysates.

Also required: Calibrator A. REF: P000057 &amp; Rinsing Fluid, REF. 8002171

## Measuring principle

Pyruvate is enzymatically oxidized by pyruvate oxidase (PyrOx). The hydrogen peroxide formed reacts with N-ethyl-N-(2-hydroxy-3-sulfopropyl)-m-toluidine (TOOS) and 4-aminoantipyrine. This reaction is catalyzed by peroxidase (POD) and yields the red-violet colored quinonediimine. The rate of formation is measured photometrically at 530 nm and is proportional to the pyruvate concentration.

Default linear range : 10 - 300 µmol/L

## Reagents:

1. Reagent: 5 bottles of lyophilisate

2. Buffer: 5 bottles of 6 mL

Reagent sufficient for 5 x 350 determinations.

Reagents are stable up to expiry date when stored at +2 to +8°C.

## Preparation and stability of solution

1. Unscrew the cap with the membrane from the reagent bottle. Remove and discard the rubber stopper.

2. Transfer the contents of the buffer bottle to the reagent bottle.

3. Fasten the cap with the membrane on the reagent bottle, without Rubber stopper.

4. Dissolve contents completely by gently turning the bottle upside-down at least ten times. Let the reagent stand and equilibrate in room temperature for at least 30 minutes prior to use.

Reconstituted reagent is stable for five days in the instrument.

	Component	Concentration in test solution
Pyruvate reagent	4-amino-antipyrine	0.3 mmol/L
	Tiamine pyrophosphate	0.2 mmol/L
	FAD	10 µmol/L
	Pyruvate Oxidase	> 0.2 kU/L
	Peroxidase	> 0.8 kU/L
Pyruvate buffer	Ascorbate Oxidase	> 10 kU/L
	Citrate buffer, pH 6.1	100 mmol/L
	Potassium dihydrogenphosphate	10 mmol/L
	MgCl <sub>2</sub>	10 mmol/L
	TOOS	1.5 mmol/L
	Sodium azide	0.3 g/L

Sample material  
Microdialysates

Calibration  
Use Calibrator A. REF. P000057

For in vitro use only

## Symbol declaration:

 Last day of use

 Lot number

 Storage temperature

 See instructions for use

 In vitro diagnostic reagent

 The product meets EU directive for IVD (98/79/EC)/LVFS 2001:7

## WARNING:

Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.

The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.

Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide

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GB: Quality Control:

It is recommended that the system be controlled using Control Samples. These controls should be assayed according to local quality assurance schemes. The use of other control materials has not been evaluated. Satisfactory level of performance is achieved when the analyte values for the controls are within the "Acceptable Control Range" published in the Packet Insert with the controls. Additional information is found in the Technical Manual for the Microdialysis Analyzers.

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M Dialysis AB