# REAGENT Set B

## For ISCUS/ISCUS<sup>flex</sup> Clinical Microdialysis Analyzer

Ref. No. 8002164

Intended user: Medical or laboratory professional staff  
Content

1. Reagent: One bottle lyophilized reagent each for Glucose, Lactate, Pyruvate, Glycerol and Glutamate placed in a cassette.
2. Buffer: One bottle 6 mL each for Glucose, Lactate, Pyruvate and Glycerol and 4 mL for Glutamate.
3. Calibrator: One bottle of Calibrator A 6 mL placed in the Reagent Cassette.

### Preparation

1. Unscrew the cap with the membrane from the reagent and calibrator bottles, placed in the cassette. Remove and discard the rubber stoppers.
2. Unscrew the cap from the buffer bottles and gently transfer the contents to the corresponding reagent bottle.
3. Fasten the cap with the membrane on the reagent and calibrator bottles in the cassette, without returning the rubber stoppers.
4. Let the reagents and calibrator stand and equilibrate in room temperature for at least 30 minutes prior to use. The contents will be completely mixed when the Reagent Cassette is placed and identified in the analyzer.

### Intended purpose and stability of solution


The Reagent Set B is a cassette containing reagents and Calibrator A made for use in ISCUS/ISCUS<sup>flex</sup> Microdialysis Analyzer. For intended purpose, see information for the individual components. Each Reagent Cassette has a unique code, written on the label, which should be registered when placing it in the analyzer.

The Reagent Set is stable up to expiry date when stored at +2 to +8 °C. Reconstituted reagents are stable for five days in the analyzer. The contents are sufficient for around 215 analyses.

## GLUCOSE

Colorimetric method for the quantitative determination of Glucose in Microdialysates.

**Measuring principle**  
Glucose is enzymatically oxidized by glucose oxidase (GOD). The hydrogen peroxide formed reacts with phenol and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields a red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the glucose concentration.

Linear range: 0.1 - 25 mmol/L

	Component	Concentration in test solution
Glucose reagent	4-aminoantipyrine	0.77 mmol/L
	Ascorbate oxidase	>3 kU/L
	Glucose oxidase	>1.5 kU/L
	Peroxidase	>1.5 kU/L
Glucose buffer	Phosphate buffer, pH 7.0	0.1 mol/L
	Phenol	11 mmol/L
	Sodium azide	0.4 g/L


#### References:

1. Barhem and P. Trinder, Analyst 97 (1972) 142

## LACTATE

Colorimetric method for the quantitative determination of Lactate in Microdialysates.

**Measuring principle**  
Lactate is enzymatically oxidized by lactate oxidase. The hydrogen peroxide formed reacts with 4-chlorophenol and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields a red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the lactate concentration.

Linear range: 0.1 - 12 mmol/L

	Component	Concentration in test solution
Lactate reagent	4-aminoantipyrine	0.4 mmol/L
	Lactate oxidase	>500 U/L
	Peroxidase	>500 U/L
Lactate buffer	Ascorbate oxidase	>12.0 kU/L
	PIPES buffer, pH 6.8	100 mmol/L
	4-Chlorophenol	5.4 mmol/L
	Sodium oxalate	7.5 mmol/L
	EDTA-disodium salt	5 mmol/L
	Sodium azide	0.3 g/L



#### References:

1. N. Shimojo et al. Clin Chem 35 (1989) 1992
2. H.F. Kühnle et al., J.Clin Chem BioChem 15 (1977) 171
3. T.O. Kleine et al., Dtsch Med Wschr 104 (1979) 553

## PYRUVATE

Colorimetric method for the quantitative determination of Pyruvate in Microdialysates.

**Measuring principle**  
Pyruvate is enzymatically oxidized by pyruvate oxidase (PyrOx). The hydrogen peroxide formed reacts with N-ethyl-N-(2-hydroxy-3-sulfo-propyl)-m-toluidine and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields a red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the pyruvate concentration.

Default linear range : 10 - 300 µmol/L

	Component	Concentration in test solution
Pyruvate reagent	4-amino-antipyrine	0.3 mmol/L
	Tiamine pyrophosphate	0.2 mmol/L
	FAD	10 µmol/L
	Pyruvate Oxidase	>0.2 kU/L
	Peroxidase	>0.8 kU/L
Pyruvate buffer	Ascorbate Oxidase	>10 kU/L
	Citrate buffer, pH 6.1	100 mmol/L
	Potassium dihydrogenphosphate	10 mmol/L
	MgCl <sub>2</sub>	10 mmol/L
	TOOS	1.5 mmol/L
	Sodium azide	0.3 g/L

#### References:



1. B. Sedewitz, et al., J. Bacteriol, 160 (1984) 273-278
2. M. Nawata, et al., Anal Biochem, 190 (1990) 84-87
3. H. Araki and M. Yamada, in: H. U. Bergmeyer (Editor), Methods of Enzymatic Analysis, 3rd ed., Vol 6, Verlag Chemie, Weinheim, 1984

## GLYCEROL

Colorimetric method for the quantitative determination of Glycerol in Microdialysates.

#### Measuring principle

Glycerol is phosphorylated by adenosin triphosphate (ATP) and glycerol kinase (GK) to glycerol-3-phosphate, which is subsequently oxidized in the presence of glycerol-3-phosphate oxidase (GPO). The hydrogen peroxide formed reacts with 3,5-dichloro-2-hydroxy-benzene sulphonic acid (DCHBS) and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields a red-violet colored quinoneimine (ACBS). The rate of formation is measured photometrically at 530 nm and is proportional to the glycerol concentration.

Linear range: 10 - 1500 µmol/L

	Component	Concentration in test solution
Glycerol reagent	4-aminoantipyrine	0.4 mmol/L
	ATP	1.0 mmol/L
	Glycerol kinase	>400 U/L
	Glycerol-3-phosphate-oxidase	>1.5 kU/L
	Peroxidase	>1 kU/L
	Ascorbate oxidase	>7.0 kU/L
Glycerol buffer	PIPES buffer, pH 7.6	40 mmol/L
	DCHBS	1.5 mmol/L
	Magnesium ions	17.5 mmol/L
	Sodium azide	0.2 g/L

#### References:



1. K.J. Foster and K.G.M.M. Alberti, Clin Chem 24 (1978) 1568

## GLUTAMATE

Colorimetric method for the quantitative determination of Glutamate in Microdialysates.

#### Measuring principle

Glutamate is enzymatically oxidized by glutamate oxidase (GluOx). The hydrogen peroxide formed reacts with N-ethyl-N-(2-hydroxy-3-sulfo-propyl)-m-toluidine (TOOS) and 4-amino-antipyrine. This reaction is catalyzed by peroxidase (POD) and yields the red-violet colored quinoneimine. The rate of formation is measured photometrically at 530 nm and is proportional to the glutamate.

Linear range: 1 - 150 µmol/L

	Component	Concentration in test solution
Glutamate reagent	4-amino-antipyrine	0.3 mmol/L
	Glutamate Oxidase	>0.25 kU/L
	Peroxidase	>0.8 kU/L
Glutamate buffer	Ascorbate Oxidase	>17.5 kU/L
	PIPES buffer, pH 6.8	100 mmol/L
	TOOS	1.5 mmol/L
	Sodium azide	0.3 g/L

#### References:

1. H. Kusakabe et al., Agric. Biol. Chem., 47 (1983) 179
2. H. Kusakabe et al., Agric. Biol. Chem., 47 (1983) 1323
3. A. Böhmer et al., Eur. J. Biochem. 182 (1989) 327-332
4. U. Wollenberger et al., Biosensors 4 (1989) 381-39

## CALIBRATOR A

#### Calibration values

Glucose	5.55 mmol/L
Lactate	2.5 mmol/L
Glycerol	475 µmol/L
Pyruvate	250 µmol/L
Glutamate	25 µmol/L

#### Sample material

Microdialysates

#### WARNING

Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.

The buffer contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.

Sodium Azide may react with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents, flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10 % sodium hydroxide.

#### Symbol declaration

Last day of use

Lot number

Storage temperature

Consult instructions for use

In vitro diagnostic device

For in vitro use only

The product meets EU directive for IVD (98/79/EC) / LVFS 2001:7

Manufactured by:  
M Dialysis AB  
Hammarby Fabriksväg 43  
SE-120 30 • Stockholm • Sweden  
Tel: +46-8-470 10 20  
Fax: +46-8-470 10 55  
E-mail: info@mdialysis.com  
www.mdialysis.com

USA office:  
M Dialysis Inc.  
73 Princeton Street  
N.Chelmsford • MA 01863 • USA  
Phone: +1 978 251 1940, +1-866-868-9236  
Fax: +1 978 251-1960  
E-mail: usa@mdialysis.com

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