
Instructions for use

GB

61 Microdialysis Catheter & 61 High Cut Off Microdialysis Catheter


INTENDED USE


61 Microdialysis Catheter and 61 High Cut Off Microdialysis Catheter are intended to enable microdialysis in hepatic tissue or in the intraperitoneal cavity.

USAGE IN HEPATIC TISSUE

The following procedure should be performed by a surgeon under aseptic conditions.


1.


2.

1. Introduce the tunneling needle lateral to the rectus muscle, preferably in the linea semilunaris, from the inside of the abdominal wall going out through the skin. Observe for bleeding. It is recommended that the tunneling needle is inserted as stated above on the left side of the patient's abdomen. This will help to prevent the catheter from kinking and facilitate removal.

2. Insert the catheter, with the protective tubing on, through the tip of the tunneling needle. Withdraw the tunneling needle.

3.

4.



3. Place a **6-0 resorbable** suture in the falciform ligament. **NOTE: Avoid using other sutures as they may compromise the safety of the catheter.**

4. Break the handles of the Splitable introducer before insertion to avoid damage of the hepatic parenchyma. Insert the introducer in the middle of the falciform ligament, in direction to segment IV and VIII (2-3 cm above the line between segment V and VIII) as close to the suture in the ligament as possible. Insert the introducer completely to the handle of the splitable tube.

5. Remove the introducer needle carefully while keeping the splitable tube in position. Do not reinsert the introducer needle after withdrawal. Observe for bleeding. **If an extensive bleeding occurs, remove the introducer.**



GB

6. Pull the protection tube of the catheter. **NOTE: Be very careful since the white membrane at the tip of the catheter is very fragile.**


7. Carefully insert the Microdialysis Catheter into the splitable tube. **Do not touch the white membrane of the catheter.**

Two people are required to follow the instructions under number 8 and 9 below:

8. Hold the Microdialysis Catheter in this position with a forceps around the tubing. Grasp the tabs of the splitable tube (three hands needed altogether) and pull the tabs apart, away from the indwelling catheter, until the tube splits down to its entire length. After this is accomplished, dispose the plastic strips at once.

9. **Keep holding the catheter with the forceps, otherwise it may come out of the liver.** Tie the suture already placed in the falciform ligament to the suture of the catheter. Tie the catheter as tight as possible to the falciform ligament. Cut off the excessive suture.

10. Place the fixation device (white plastic wing) close to the insertion site and secure it to the catheter by suturing over the groves/indentations. Suture the fixation device to the skin through the two holes.


11. Put a protective bandage over the insertion site, e.g. Tegaderm®, OpSite® or similar.

12. After closing of the abdomen connect the inlet tubing of the catheter to the 106 Syringe, filled with Perfusion Fluid. Place the syringe in the 106 or 107 Microdialysis Pump and close the lid. Insert a Microvial into the microvial holder.


USAGE IN INTRAPERITONEAL CAVITY


The following procedure should be performed by a surgeon under aseptic conditions. Before the abdominal wall is closed:

- 1a** **1b**
1. a,b Introduce the tunneling needle lateral to the rectus muscle, in the linea semilunaris, from the inside of the abdominal wall going out through the skin. Observe for bleeding.

- 2** **3**
2. Insert the catheter, with the protective tubing on, through the tip of the tunneling needle.
 3. Withdraw the tunneling needle.

- 4** **6**
4. Remove the protective tubing by holding the liquid cross. Unscrew the protective tubing by turning it counter clockwise.
 - NOTE:** Be sure to straighten the protective tubing before removing it from the catheter. Use caution when removing the protective tubing so that you do not damage the membrane at the distal end of the catheter.
 5. Place the membrane of the catheter close to the anastomosis or in the region of interest. The catheter shall be placed free floating in the intraperitoneal cavity.
 6. Place the fixation device over the liquid cross and secure it to the catheter by suturing around the groves/indentations on the device.
 7. Suture the fixation device to the patient's skin through the two holes located on both sides of the fixation device.
 8. Fix the tubing to the skin with a protective bandage, e.g. Tegaderm®, Opsite® or similar.
 9. Insert a microvial into the microvial holder.
 10. Connect the inlet tubing to the syringe, filled with Perfusion Fluid.
 11. Place the syringe in the 106 or 107 Microdialysis Pump and close the lid.

REMOVAL OF CATHETER

The catheter is removed by gently pulling it out through the insertion site. Make sure that the suture is still attached to the catheter.

NOTE: Discontinued/removed catheters shall be handled according to the hospital routines for biohazard material.

INDICATIONS

To monitor metabolism intraperitoneally or in hepatic tissue. When used intraperitoneally: the catheter should be placed freely in the intraperitoneal cavity, it must not be inserted or attached to the intestinal wall.

Microdialysis shall not be used as the sole basis for diagnosis or therapy.

CONTRAINDICATIONS & RISKS


- Inserting the catheter may cause bleeding from damaged vessels.
- Usage in Hepatic tissue: Do not use the catheter if the mass of the liver tissue is very small
- The 61 Microdialysis Catheter & 61 High Cut Off Microdialysis Catheter is for single use only. If the device is re-used there is a risk for cross-contamination.
- For 61 High Cut Off Microdialysis Catheter: Patients with known hypersensitivity to Dextran.

GB

PRECAUTIONS


- This device is sterile unless the blister package has been opened or damaged.
- The 61 Microdialysis Catheter & 61 High Cut Off Microdialysis Catheter shall only be used together with the accessories described in the table below.
- Be sure to handle the catheter carefully to avoid kinking or other damage, particularly after removal of the protection tube. **Avoid contact with the dialysis membrane.**
- If any visible damage is observed the catheter shall not be used.
- If there is a suspicion that the catheter has become unsterile prior to insertion the catheter shall not be used.
- The pump syringe connected to the catheter should not manually be flushed since that could damage the dialysis membrane
- Check that liquid is being pumped through the catheter by inspecting the volume in the microvial **each time** the microvials are changed.
- If there is no fluid in the collected vial, start a flush on the pump: Open the lid, wait 3 seconds and close it again. Wait for the flush (5 minutes). Check that the tubing's are not kinked and that the microvial holder needle is correctly piercing the microvial membrane. If there is still no fluid in the collected vial, the dialysis membrane might be damaged, and the catheter has to be removed.
- Remove the catheter if there is a permanent stop in the liquid flow.
- 61 Microdialysis Catheter & 61 High Cut Off Microdialysis Catheter is biocompatible up to 30 days. It may though stop working earlier because of clogging, duration of use is up to 10 days referring to literature.
- The splittable Introducer is only for use when inserting the catheter into hepatic tissue.

PRODUCT DESCRIPTION & TECHNICAL INFORMATION

61 Microdialysis Catheter & 61 High Cut Off Microdialysis Catheter	REF 8010226	REF 8050191
1. Microvial (Polystyrene + Santoprene)	-	-
2. Luer lock connection (Polycarbonate)	-	-
3. Inlet tube (Polyurethane, OD 1 mm)	400 mm	400 mm
4. Vial holder (Polycarbonate)	-	-
5. Outlet tube (Polyurethane, OD 1 mm)	70 mm	70 mm
6. Liquid cross (Polysulfone)	-	-
7. Protection tube (Polyethylene)	-	-
8. Inner Shaft (Polyurethane, OD 0.9 mm)	310 mm	310 mm
9. Outer Shaft (Polyurethane, OD 1.5 mm)	280 mm	280 mm
10. Suture (Braided polyester suture, 5-0, non-resorb able)	approx. 200 mm	approx. 200 mm
11a. Dialysis membrane (Polyarylethersulphone, PAES, OD 0.6 mm)	30 mm	30 mm
11b. Membrane Cut off	20 000 Dalton	100 000 Dalton
12. Gold thread within the catheter membrane tip (OD 0.13 mm)	3 mm	3 mm
13. Fixating device (PEBA X®)	-	-
Dead Volume Time - From membrane to microvial (at 0,3 µl/min)	24 min	24 min

The distal part of the catheter has a gold thread within the catheter tip, which makes it possible to detect the location of the catheter in the tissue through CT-scanning.

GB

SPLITABLE INTRODUCER SI-2	REF 8010343
14. Handle, Cannula (Polycarbonate)	-
15. Handle, Splitable tube (Polytetrafluorethylene)	-
16. Splitable tube, (Polytetrafluorethylene, OD 1.5 mm)	74 mm
17. Cannula (Stainless steel, OD 1.0 mm)	77 mm

ACCESSORIES

The 61 Microdialysis Catheter & 61 High Cut Off Microdialysis Catheter shall only be used with these accessories.

REF	Name
P000034	Perfusion fluid T1 10x5 mL
8010191	106 Syringe 20/pkg
P000154	Microvials in rack, Sterile 12x4
P000055	Tunneling needle
P000003	106 Microdialysis Pump
P000127	107 Microdialysis Pump

CONSUMABLES

The 61 Microdialysis Catheter & 61 High Cut Off Microdialysis Catheter should be used with these consumables.

REF	Name
P000001	Microvials 250/pkg

For 61 High Cut Off Microdialysis Catheter:

To avoid ultra-filtration over the membrane the catheter should be used with a Perfusion Fluid consisting of a Ringer solution containing 30g Dextran 60/1000mL or similar.

SYMBOL EXPLANATION & PRODUCT LIMITATIONS


	Last date of use (YYYY-MM-DD)
	LOT number
	Catalogue number
	Single use only
	See instructions for use
	Fulfills the Medical Device Directive 93/42/EEC.
	Sterilised by β -radiation
	Storage temperature (4-25° C)
	Manufacturer
	Do not use if package is damaged
Membrane	Membrane material (Polyarylethersulphone, PAES)
Memb. length	Membrane length (mm)
Shaft length	Shaft length (mm)
Outlet	Outlet length (mm)
Inlet	Inlet length (mm)

Bruksanvisning

SE

61 Microdialysis Catheter & 61 High Cut Off Microdialysis Catheter

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