
Instructions for use

GB

62 GASTROINTESTINAL MICRODIALYSIS CATHETER

INTENDED USE

The 62 Gastrointestinal Microdialysis Catheter is a single use product that is intended to measure local tissue chemistry by performing microdialysis on patients undergoing gastrointestinal surgery.

PRODUCT DESCRIPTION & TECHNICAL INFORMATION

| | REF 8010292 |
|--|------------------------|
| 1. Microvial (Polystyrene + Santoprene) | X |
| 2. Luer lock connection (Polycarbonate) | X |
| 3. Inlet tube (Polyurethane, OD 1mm) | 600 mm |
| 4. Vial holder (Polycarbonate) | X |
| 5. Outlet tube (Polyurethane, OD 1mm) | 220 mm |
| 6. Fixation device (PEBAX®) | X |
| 7. Liquid cross (Polysulfone) | X |
| 8. Protection tube (Polyethylene) | 240 mm |
| 9. Shaft (Polyurethane, OD 0.9 mm) | 180 mm |
| 10. Dialysis membrane (Polyarylethersulfone, PAES OD 0.6 mm) Cut off:20 000 Dalton | 30 mm |
| 11. Gold thread within the catheter membrane tip (OD 0.13 mm; L=3 mm) | X |
| Dead Volume Time - From membrane to microvial (at 0,3 µl/min) | 24min |

The distal part of the catheter has a gold thread within the catheter tip, which makes it possible to detect the location of the catheter in the tissue through CT-scanning.

ACCESSORIES

The 62 Gastrointestinal Microdialysis Catheter shall only be used with these accessories.

| REF | Name |
|------------|----------------------------------|
| P000034 | Perfusion fluid T1 10x5 mL |
| 8010191 | 106 Syringe 20/pkg |
| P000154 | Microvials in rack, Sterile 12x4 |
| P000055 | Tunneling needle |
| P000003 | 106 Microdialysis Pump |
| P000127 | 107 Microdialysis Pump |

CONSUMABLES

The 62 Gastrointestinal Microdialysis Catheter should be used with these consumables.

| REF | Name |
|------------|--------------------|
| P000001 | Microvials 250/pkg |

INDICATIONS

The 62 Gastrointestinal Microdialysis Catheter is designed to be placed free floating in the intraperitoneal cavity. It shall not be introduced or attached into the intestinal wall. Microdialysis shall not be used as the sole basis for diagnosis or therapy.

CONTRAINDICATIONS & RISKS

- Inserting the catheter may cause bleeding from damaged vessels.
- The 62 Gastrointestinal Microdialysis Catheter is for single use only. If the device is re-used there is a risk for cross-contamination.



PRECAUTIONS



- This device is sterile unless the package has been opened or damaged.
- The 62 Gastrointestinal Microdialysis Catheter shall only be used together with the accessories described in the previous table.
- Be sure to handle the catheter carefully to avoid kinking or other damage, particularly after removal of the protection tube. Avoid contact with the dialysis membrane.
- If any visible damage is observed the catheter shall not be used.
- If there is a suspicion that the catheter has become unsterile prior to insertion the catheter shall not be used.
- The pump syringe connected to the catheter should not manually be flushed since that could damage the dialysis membrane.
- Check that liquid is being pumped through the catheter by inspecting the volume in the microvial **each time** the microvials are changed.
- If there is no fluid in the collected vial, start a flush on the pump: Open the lid, wait 3 seconds and close it again. Wait for the flush (5 minutes). Check that the tubing's are not kinked and that the microvial holder needle is correctly piercing the microvial membrane. If there is still no fluid in the collected vial, the dialysis membrane might be damaged and the catheter has to be removed.
- Remove the catheter if there is a permanent stop in the liquid flow.
- 62 Gastrointestinal Microdialysis Catheter is biocompatible up to 30 days. It may though stop working earlier because of clogging, duration of use is up to 10 days referring to literature.

USAGE



The following procedure should be performed by a surgeon under aseptic conditions. Before the abdominal wall is closed:

- 1a, b** Introduce the tunneling needle lateral to the rectus muscle, in the linea semilunaris, from the inside of the abdominal wall going out through the skin. Observe for bleeding.

- 2** Insert the catheter, with the protective tubing on, through the tip of the tunneling needle.
3 Withdraw the tunneling needle.

- 4** Remove the protective tubing by holding the liquid cross. Unscrew the protective tubing by turning it counter clockwise.
NOTE: Be sure to straighten the protective tubing before removing it from the catheter. Use caution when removing the protective tubing so that you do not damage the membrane at the distal end of the catheter.
5 Place the membrane of the catheter close to the anastomosis or in the region of interest. The catheter shall be placed free floating in the intraperitoneal cavity.
6 Place the fixation device over the liquid cross and secure it to the catheter by suturing around the groves/indentations on the device.
7 Suture the fixation device to the patient's skin through the two holes located on both sides of the fixation device.
8 Fix the tubing to the skin with a protective bandage, e.g. Tegaderm®, Opsite® or similar.
9 Insert a microvial into the microvial holder.
10 Connect the inlet tubing to the syringe, filled with Perfusion Fluid T1.
11 Place the syringe in the 106 or 107 Microdialysis Pump and close the lid.

REMOVAL OF CATHETER

The catheter is removed by gently pulling it out through the insertion site.

GB

NOTE: Discontinued/removed catheters shall be handled according to the hospital routines for biohazard material.

SYMBOL EXPLANATION & PRODUCT LIMITATIONS

| | |
|--------------|--|
| | Last date of use (YYYY-MM-DD) |
| | LOT number |
| | Catalogue number |
| | Single use only |
| | See instructions for use |
| | Fulfils the Medical Device Directive 93/42/EEC |
| | Sterilised by β -radiation |
| | Storage temperature (4-25° C) |
| | Manufacturer |
| | Do no use if package is damaged |
| Membrane | Membrane material (Polyarylethersulfone, PAES) |
| Memb. length | Membrane length (mm) |
| Shaft length | Shaft length (mm) |
| Outlet | Outlet length (mm) |
| Inlet | Inlet length (mm) |

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