
## Instructions for use

### 66 Linear Catheter & 66 High Cut Off Linear Catheter

*For precautions and further details, see next page.*


**INTENDED USE**

The **66 Linear Catheter** and **66 High Cut Off Linear Catheter** is intended to enable microdialysis in skin, adipose and resting skeletal muscle tissues.


**USAGE**

The procedure for insertion shall be performed by a physician or a nurse under aseptic conditions. The following procedure describes insertion in skin tissue.

1. Make marks for catheter entrance and exit holes with circles (Fig 1).
2. Take the introducer cannula and insert it through the circle marking the catheter exit hole (Fig 2) and continue through the skin to the marked catheter entrance hole (Fig 3). Avoid penetration of actual markings.


3. Remove the catheter from the protective tube and guide the catheter gently in to the sharp end of the introducer cannula (Fig 4). Push it **carefully** through the introducer cannula until it appears at the opposite side of the cannula and continue until the membrane is in position (Fig 5).  
**Important:** Be very gentle when entering the membrane into the introducer cannula.
4. When the catheter and membrane is in position, the introducer cannula is removed (Fig 6), after which the catheter is fixed to the skin.


5. Connect the sample vial (Fig 7, 8, 9). It is important to position the vial at the same level or lower than the membrane in order to avoid ultrafiltration through the membrane.


**Fig 7 (a, b)**

**Fig 8 (a, b, c)**

**Fig 9 (a, b)**

6. Connect the catheter to the syringe filled with perfusion fluid (106 Syringe) and place the syringe in the pump (106 or 107) and close the lid. Check that fluid comes out of the outlet, and into the vial within approximately 2 minutes.

**REMOVAL OF THE CATHETER**

The removal of the catheter shall be performed under aseptic conditions. Cut the outlet tube as short as possible. The catheter is then removed by carefully pulling it out of the tissue.  
**NOTE:** Discontinued/removed catheters shall be handled according to the hospital routines for biohazard material.

**CONTRAINDICATIONS & RISKS**

For **66 High Cut Off Linear Catheter**: Patients with known hypersensitivity to Dextran. The 66 Linear Catheter and 66 High Cut Off Linear is for single use only. If the device is re-used there is a risk for cross-contamination. Microdialysis should not be used as the sole basis for diagnosis or therapy.

**PRECAUTIONS**

- Care must be taken to avoid needle stick injuries to reduce risk of exposure to contaminated blood.
- The 66 Linear Catheter and 66 High Cut Off Linear shall only be used with the introducer cannula included in the package.
- To avoid ultra-filtration over the membrane the 66 **High Cut Off Linear Catheter** should be used with a Perfusion Fluid consisting of a Ringer solution containing 30g Dextran 60/1000mL or similar.
- Remove the catheter if there are any signs of infection.
- Utstyret er sterilt så lenge blisterpakken ikke har vært åpnet eller er skadet.
- 66 Linear Catheter og 66 High Cut Off Linear skal kun brukes sammen med tilbehøret som er beskrevet i tabellen nedenfor.
- Be sure to handle the catheter carefully to avoid kinking or other damage, particularly after removal of the protection tube. **Avoid contact with the dialysis membrane.**
- If any visible damage is observed, do not use the catheter.
- If there is a suspicion that the catheter has become unsterile prior to insertion the catheter shall not be used.
- The pump syringe connected to the catheter should not manually be flushed since that could damage the dialysis membrane.
- Check that liquid is being pumped through the catheter by inspecting the volume in the microvial **each time** the microvials are changed.
- If there is no fluid in the collected vial, start a flush on the pump: Open the lid, wait 3 seconds and close it again. Wait for the flush (5 minutes). Check that the tubing's are not kinked. If there is still no fluid in the collected vial, the dialysis membrane might be damaged, and the catheter has to be removed.
- Remove the catheter if there is a permanent stop in the liquid flow.


**DELAY IN THE SYSTEM AND THE TIME IT TAKES TO FILL UP A MICROVIAL**

Flow rate	Delay time (with 100 mm outlet)	Time to acquire minimum volume for analysis in analyzers	Time to fill a microvial (200 μl)
0,3 μl/min	≈ 6 min	9 – 14 min*	10 hours
0,5 μl/min	≈ 3 min	5 – 8 min*	6 hours
1 μl/min	≈ 2 min	3 – 4 min*	3 hours

\*Depending on the number of substances that shall be analyzed

**PRODUCT DESCRIPTION & TECHNICAL INFORMATION**

X= Identical on all catheters

	REF 8010650	REF 8010651 (High cut off)	REF 8010670	REF 8010671 (High cut off)
1. Microvial with vial cap (Polystyrene + Santoprene)	X	X	X	X
2. Luerlock-connector (Polycarbonate)	X	X	X	X
3. Inlet tube (Polyurethane, OD 1 mm)	50 mm	50 mm	50 mm	50 mm
4. Inlet tube (Polyurethane, OD 0,4 mm)	350 mm	350 mm	350 mm	350 mm
5. Membrane (Polyarylethersulphone, PAES, OD 0.5 mm)	30 mm	30 mm	10 mm	10 mm
6. Outlet tube (Polyurethane, OD 0,4 mm)	100 mm	100 mm	100 mm	100 mm
7. Protective tube (Polyethylene)	145 mm	145 mm	145 mm	145 mm
8. Introducer cannula 21G (Stainless steel))	50 mm	50 mm	50 mm	50 mm
9. Vial Cap	X	X	X	X
Membrane Cut off (Dalton)	20 000	100 000	20 000	100 000

**ACCESSORIES**

66 shall only be used with the following accessories:

REF	Name
8010191	106 Syringe 20/pkg
P000154	Microvials in rack, Sterile 12x4
P000003	106 Microdialysis Pump
P000127	107 Microdialysis Pump
P000034	Perfusion Fluid T1

**CONSUMABLES**

The 66 should be used with these consumables.

REF	Name
P000001	Microvials 250/pkg
7431007	CMA Vial Glass 300 μl (supplied by CMA Microdialysis AB)

**For 66 High Cut Off Linear Catheter:**

To avoid ultra-filtration over the membrane the catheter should be used with a Perfusion Fluid consisting of a Ringer solution containing 30g Dextran 60/1000mL or similar.


**GB**

**SYMBOL EXPLANATION & PRODUCT LIMITATIONS**

	Last month of use (YYYY-MM-DD)
	LOT number
	Catalogue number
	Single use only
	See instructions for use
	Fulfils the Medical Device Directive 93/42/EEC.
	Sterilised by β-radiation
	Storage temperature (4-25° C)
	Manufacturer
	Do not use if package is damaged
Membrane	Membrane material, (Polyarylethersulphone, PAES)
Memb. length	Membrane length (mm)
Outlet	Outlet length (mm)
Inlet	Inlet length (mm)

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