



- Perfusion fluid for microdialysis in the brain
- Easy and safe to use

## Perfusion Fluid CNS for Microdialysis

Perfusion Fluid CNS is an isotonic sterile perfusion fluid especially developed for Microdialysis in CNS tissue. It is conveniently packed in glass ampoules, 5 mL, ready to use. Perfusion Fluid CNS is cleared as a medical device according to the CE medical directive.

### Ordering information

**Features** Isotonic, sterile  
**Contents** NaCl 147 mmol/L  
 KCl 2.7mmol/L  
 CaCl<sub>2</sub> 1.2 mmol/L  
 MgCl<sub>2</sub> 0.85 mmol/L

---

Total Chloride contents:  
 153.8 mmol/L

**Osmolality** 305 mOsm/kg  
**pH** ~6

**Package** Glass ampoules, 10 x 5 mL  
**Ref. No.** P000151, Perfusion fluid CNS

**CE** CE marked according to the Medical Device Directive,93/42/EEC

**STERILE** Sterilized by steam

**2** Single use only

 **dialysis**