



- Perfusion fluid for microdialysis
- Easy and safe to use

Perfusion Fluid T1 for Microdialysis

Perfusion Fluid T1 is an isotonic sterile perfusion fluid especially developed for Microdialysis. It is conveniently packed in glass ampoules, 5 mL, ready to use. Perfusion Fluid T1 is cleared as a medical device according to the CE directive.

Ordering information

Features Isotonic, sterile
Contents NaCl 147 mmol/L
 KCl 4 mmol/L
 CaCl₂ 2.3 mmol/L

Total Chloride contents:
 155.6 mmol/L

Osmolality 290 mOsm/kg
pH ~6

Package Glass ampoules, 10 x 5 mL
Ref. No. P000034, Perfusion fluid T1



CE marked according to the Medical Device Directive, 93/42/EEC



Sterilized by steam



Single use only

