



- **Perfusion fluid for use with 71 High Cut off Brain Microdialysis catheter**

Perfusion Fluid CNS Dextran for Microdialysis

Perfusion Fluid CNS Dextran is an isotonic sterile perfusion fluid especially developed for Microdialysis in brain tissue with the 71 High Cut off Brain Microdialysis Catheter. It is conveniently packed in glass bottles, 5 mL, ready to use. Perfusion Fluid CNS Dextran is cleared as a medical device according to the CE medical directive.

Ordering information

Features: Isotonic, sterile

Contents: NaCl 147 mmol/L

KCl 2.7mmol/L

CaCl₂ 1.2 mmol/L

MgCl₂ 0.85 mmol/L

Total Chloride contents: 153.8 mmol/L

Dextran MW 500 kDa 3%

pH: ~6

Package Glass bottles, 2 x 5 mL

Ref. No. 8050151, Perfusion fluid CNS Dextran



CE marked according to the Medical Device Directive, 93/42/EEC



Sterilized by use of aseptic processing technique



Single use only

 **dialysis**