

# CONTROL SAMPLES

## for Microdialysis Analyzers

**USA**

| LOT                              | T26842 | Component | Acceptable control range |        |
|----------------------------------|--------|-----------|--------------------------|--------|
| Control Low<br>Ref:8010403       |        | Glucose   | 19.91 - 26.93            | mg/dL  |
|                                  |        | Lactate   | 1.20 - 2.00              | mmol/L |
|                                  |        | Pyruvate  | 13.75 - 22.92            | μmol/L |
|                                  |        | Glycerol  | 48.75 - 81.25            | μmol/L |
|                                  |        | Glutamate | 7.50 - 12.50             | μmol/L |
| Control Elevated<br>Ref: 8010205 |        | Glucose   | 238.89 - 323.21          | mg/dL  |
|                                  |        | Lactate   | 7.68 - 11.52             | mmol/L |
|                                  |        | Pyruvate  | 187.00 - 253.00          | μmol/L |
|                                  |        | Glycerol  | 663.00 - 897.00          | μmol/L |
|                                  |        | Glutamate | 102.00 - 138.00          | μmol/L |

### QR Codes

**Nominal Control Levels**

Service Code      Levels



**CLIA Intervals**

Service Code



Intervals



<http://www.mdialysis.com/products/clinical-microdialysis-international/products/analyzers/controlsamplesdownload>

