

# CONTROL SAMPLES

## for Microdialysis Analyzers

**USA**

| LOT                              | T26977 | Component | Acceptable control range |
|----------------------------------|--------|-----------|--------------------------|
| Control Low<br>Ref:8010403       |        | Glucose   | 19.91 - 26.93 mg/dL      |
|                                  |        | Lactate   | 1.20 - 2.00 mmol/L       |
|                                  |        | Pyruvate  | 13.75 - 22.92 μmol/L     |
|                                  |        | Glycerol  | 48.75 - 81.25 μmol/L     |
|                                  |        | Glutamate | 7.50 - 12.50 μmol/L      |
| Control Elevated<br>Ref: 8010205 |        | Glucose   | 238.89 - 323.21 mg/dL    |
|                                  |        | Lactate   | 7.68 - 11.52 mmol/L      |
|                                  |        | Pyruvate  | 187.00 - 253.00 μmol/L   |
|                                  |        | Glycerol  | 663.00 - 897.00 μmol/L   |
|                                  |        | Glutamate | 102.00 - 138.00 μmol/L   |

### QR Codes

**Nominal Control Levels**

Service Code

Levels



**CLIA Intervals**

Service Code

Intervals



<http://www.mdialysis.com/products/clinical-microdialysis-international/products/analyzers/controlsamplesdownload>

