



- Perfusion fluid for microdialysis in the brain
- Easy and safe to use

Perfusion Fluid CNS for Microdialysis

Perfusion Fluid CNS is an isotonic sterile perfusion fluid especially developed for Microdialysis in CNS tissue. It is conveniently packed in glass ampoules, 5 mL, ready to use. Perfusion Fluid CNS is cleared as a medical device according to the CE medical directive.

Ordering information

Features Isotonic, sterile
Contents NaCl 147 mmol/L
 KCl 2.7 mmol/L
 CaCl₂ 1.2 mmol/L
 MgCl₂ 0.85 mmol/L

Total Chloride contents:
 153.8 mmol/L

Osmolality 305 mOsm/kg
pH ~6

Package Glass ampoules, 10 x 5 mL
Ref. No. P000151, Perfusion fluid CNS



CE marked according to the Medical Device Directive, 93/42/EEC



Sterilized by steam



Single use only

μ dialysis