



- Perfusion fluid for microdialysis
- Easy and safe to use

## Perfusion Fluid T1 for Microdialysis

Perfusion Fluid T1 is an isotonic sterile perfusion fluid especially developed for Microdialysis. It is conveniently packed in glass ampoules, 5 mL, ready to use. Perfusion Fluid T1 is cleared as a medical device according to the CE directive.

### Ordering information

**Features** Isotonic, sterile  
**Contents** NaCl 147 mmol/L  
 KCl 4 mmol/L  
 CaCl<sub>2</sub> 2.3 mmol/L

Total Chloride contents:  
 155.6 mmol/L

**Osmolality** 290 mOsm/kg  
**pH** ~6

**Package** Glass ampoules, 10 x 5 mL  
**Ref. No.** P000034, Perfusion fluid T1



CE marked according to the Medical Device Directive, 93/42/EEC



Sterilized by steam



Single use only

*μ* dialysis