

Magnus Hedberg 2022-03-14

ISCUS^{flex} New Software Version March 2022

- Rev L Version 2.1.0.494



ISCUS^{flex} New Software Version March 2022

- Rev L Version 2.1.0.494



Overview

- Lactate calibrator response fix
- Difficulties to measure Glucose blank. Fixed, notice in service log.
- Default sample and reagent volumes according to service code SAMPLEVOL OPTION

Defect fixes

(TT #335 – TT #338)

- TT #335 Lactate calibrator response too high. Automatic update at startup. Implemented also as an ISCUS tools fix.
- TT #336 Difficult to measure Glucose blank. Fixed, notice in service log. Handles initially noisy reagents.
- TT #337 Default sample and reagent volumes set to SAMPLEVOL OPTION
- TT #338 Faulty sample volumes for preclinical Glucose. Corrected

Service codes

- [D6558-17 P192 ISCUSflex Service codes \(Rev L 485 and later\)](#)