



◦ Perfusion fluid for use with
71 High Cut off Microdialysis catheters

Perfusion Fluid CNS Dextran for Microdialysis

Perfusion Fluid CNS Dextran is an isotonic sterile perfusion fluid especially developed for Microdialysis in brain tissue with the 71 High Cut off Brain Microdialysis Catheter and 71 High Cut Off Microdialysis Bolt Catheter. It is conveniently packed in glass bottles, 7,5 mL, ready to use.

Perfusion Fluid CNS Dextran is cleared as a medical device according to the CE medical directive.

Ordering information

Features: Isotonic, sterile

Package: Glass bottles, 2 x 7,5 mL

Contents: NaCl 147 mmol/L

Ref. No.: 8050151, Perfusion fluid CNS Dextran

KCl 2.7 mmol/L

CaCl₂ 1.2 mmol/L

MgCl₂ 0.85 mmol/L

Total Chloride contents: 153.8 mmol/L

Dextran MW 500 kDa 3%

pH:

~6

CE 0413

CE marked according to the Medical Device Directive, 93/42/EEC

STERILE A

Sterilized by use of aseptic processing technique



Single use only



Shelf life: 2 years from the manufacturing date



Storage temperature: 4-25 °C



Mikrosol AB Hammarby Fabriksväg 43, SE-120 30 Stockholm, Sweden Tel: +46-8-470 10 20 Fax: +46-8-470 10 55 E-mail: info@mdialysis.com www.mdialysis.com

US Office: 73 Princeton Street, North Chelmsford, MA 01863 USA Tel: (800) 440-9236 (978) 251-1940 Fax: (978) 251-1960 E-mail: usa@mdialysis.com

μ dialysis