



- Perfusion fluid for microdialysis
- Easy and safe to use

Perfusion Fluid T1 for Microdialysis

Perfusion Fluid T1 is an isotonic sterile perfusion fluid especially developed for Microdialysis. It is conveniently packed in glass ampoules, 5 mL, ready to use.

Intended use: The Perfusion Fluid T1 is intended to act as a carrier inside the M Dialysis Microdialysis catheters.

Ordering information


Features Isotonic, sterile
Contents NaCl 147 mmol/L
 KCl 4 mmol/L
 CaCl₂ 2.3 mmol/L


Package Glass ampoules, 10 x 5 mL
Ref. No. P000034, Perfusion fluid T1

Total Chloride contents:
 155.6 mmol/L

CE 0413 CE marked according to the Medical Device Directive, 93/42/EEC

Osmolality 290 mOsm/kg
pH -6

STERILE  Sterilized by steam

 Single use only



80038641, 2022-12-21