



- Perfusion fluid for microdialysis in the brain
- Easy and safe to use

## Perfusion Fluid CNS for Microdialysis

Perfusion Fluid CNS is an isotonic sterile perfusion fluid especially developed for Microdialysis in CNS tissue. It is conveniently packed in glass ampoules, 5 mL, ready to use.

Intended use: The Perfusion Fluid CNS is intended to act as a carrier inside the M Dialysis Microdialysis catheters.

### Ordering information


**Features** Isotonic, sterile  
**Contents** NaCl 147 mmol/L  
 KCl 2.7 mmol/L  
 CaCl<sub>2</sub> 1.2 mmol/L  
 MgCl<sub>2</sub> 0.85 mmol/L

Total Chloride contents:  
 153.8 mmol/L

**Osmolality** 305 mOsm/kg  
**pH** ~6

**Package** Glass ampoules, 10 x 5 mL  
**Ref. No.** P000151, Perfusion fluid CNS

CE 0413 CE marked according to the Medical Device Directive, 93/42/EEC

STERILE  Sterilized by steam

 Single use only

