




- Perfusion fluid for microdialysis in brain tissue
- Easy and safe to use

Perfusion Fluid CNS for Microdialysis

Perfusion Fluid CNS is an isotonic sterile perfusion fluid especially developed for Microdialysis monitoring in brain tissue by the Cerebral Tissue Monitoring system from MDialysis. The fluid is conveniently packed in glass bottles 7.5 mL, ready to use.

Ordering information

Features Isotonic, sterile
Contents NaCl 147 mmol/L
KCl 2.7 mmol/L
CaCl₂ 1.2 mmol/L
MgCl₂ 0.85 mmol/L

Total Chloride contents:
153.8 mmol/L

Osmolality 305 mOsm/kg
pH ~6

Package Glass bottles, 10 x 7.5 mL
Ref. No. P000151, Perfusion fluid CNS



Sterilized using steam



Single use only

 **dialysis**